

PRODUCT CATALOG SUPPORTS AND ORTHOSES

2020

BAUERFEIND.COM

BAUERFEIND – BY PEOPLE FOR PEOPLE. FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Our specialists have worked very closely with physicians, scientists and developers right from the start. The result: patented, award-winning products and services.

From the thread to the completed product, Bauerfeind supports and orthoses are developed and manufactured in Germany.

EFFECTIVE VARIETY

Bauerfeind has taken on the mission of using modern orthopedic supports and orthoses to help people lead an active and pain-free life, even well into old age. Our products are characterized by clinically proven effectiveness and unique material properties.

The supports and orthoses help during recovery from injuries, in cases of degenerative changes and following surgery, and they provide new levels of mobility. Bauerfeind's high quality standards guarantee excellent product functionality.

QUALITY WITHOUT COMPROMISE

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements in terms of quality, effectiveness, and wearing comfort. This is confirmed by regular inspection procedures, for which our products are permitted to bear the seals and marks after successful testing. All of Bauerfeind's supports and orthoses conform to the requirements of the German Medical Devices Act and the EU Medical Devices Directive 93/42/EEC. The CE certification and DIN EN ISO 13485 certification guarantee quality you can trust.

Please note:

Medical aids from Bauerfeind can be distributed irrespective of age, as long as the dimensions of the relevant body part fit the product's specific sizing system.









CONTENTS

THE BAUERFEIND SERVICE

YOU CAN RELY ON IT

> Page 6

BAUERFEIND PRODUCT FEATURES A PLUS IN TERMS OF FUNCTION, QUALITY AND WEARING COMFORT

>Page 8

THE NEW TRAIN SUPPORTS

THE BENCHMARK FOR THE FUTURE

> Page 10

THE BAUERFEIND MEDICAL AID APP

THE PRODUCT FINDER FOR PHYSICIANS AND MEDICAL RETAILERS

> Page 12

THE BAUERFEIND TREATMENT APP

TARGETED TRAINING FOR YOUR CUSTOMERS

> Page 13

THE RIGHT PRODUCT FOR EVERY TREATMENT AREA



FOOT AND ANKLE

- 16 Table of indications
- **18** MalleoTrain®
- **19** MalleoTrain[®] S
- **20** MalleoTrain[®] S open heel
- 21 MalleoTrain[®] Plus
- 22 AchilloTrain®
- 23 AchilloTrain® Pro
- **24** AirLoc[®]
- 25 MalleoLoc[®] L
- 26 MalleoLoc® L3
- 28 MalleoLoc®
- 29 CaligaLoc®
- **30** ValguLoc[®]
- 31 ValguLoc® II

ADDITIONAL PRODUCT INFORMATION AND SERVICES

ACCESSORIES AND SPARE PARTS

- 106 Foot and ankle
- **107** Knee, thigh and hip
- **110** Spine and neck
- 114 Hand, elbow and shoulder



KNEE, THIGH AND HIP

- 34 Table of indications
- 36 GenuTrain®
- **38** GenuTrain[®] A3
- **39** GenuTrain[®] P3
- 40 GenuTrain® S
- 41 GenuTrain[®] S Pro
- 42 GenuPoint®
- 43 SofTec® Genu
- 44 SecuTec® Genu
- 46 GenuTrain® OA
- 48 SofTec® OA
- 49 SecuTec® OA
- 50 MOS-Genu
- **51** GenuLoc[®]
- **52** CoxaTrain[®]
- 54 MyoTrain®
- 55 SofTec® Coxa

50

SPINE AND NECK

58 Table of indications

- 60 LumboTrain[®] | Lady
- 62 SacroLoc®
- 64 LordoLoc®
- 65 LumboLoc[®]
- 66 LumboLoc® Forte
- **67** Spinova[®] system
- 70 Spinova® Stabi Classic
- 71 Spinova[®] Support Plus
- 72 Spinova® Support Plus Classic
- 73 Spinova® Unload Classic
- 74 Spinova[®] Unload Plus
- 75 Spinova[®] Immo Plus
- 76 Spinova® Immo Plus Classic
- 77 Spinova[®] Immo
- 78 SecuTec® Dorso
- 79 SecuTec® Lumbo
- 80 SofTec® Lumbo
- 81 DorsoTrain®
- 82 Spinova® Osteo
- 84 SofTec® Dorso
- 84 CerviLoc[®] | CerviLoc[®] S

HAND, ELBOW AND SHOULDER

- 88 Table of indications
- 90 ManuTrain®
- 91 ManuLoc®
- 92 ManuLoc® long
- **93** ManuLoc[®] long Plus
- 94 RhizoLoc®
- 95 ManuLoc® Rhizo
- 96 ManuLoc® Rhizo long
- 97 ManuLoc[®] Rhizo long Plus
- 98 EpiTrain®
- 99 EpiPoint®
- 100 OmoTrain®
- 101 OmoTrain[®] S
- 102 OmoLoc®
- 103 SecuTec® Omo

CUSTOM-MADE AND STANDARD PRODUCTS

- **118** Order sheets for custom-made products
- **125** Order sheet for standard products
- **126** Bauerfeind measurement technology

THE BAUERFEIND SERVICE YOU CAN RELY ON IT

ONLINE SHOP SERVICE

When ordering via the Bauerfeind online shop, you benefit from reduced shipping costs and can take advantage of versatile functions such as the order status, availability display or patient files. To create a customer account, register at: **experts.bauerfeind.com/register**



YOUR BENEFITS AT A GLANCE

- + 24 h online availability
- + Accessible anywhere using a web browser
- + Intuitive order process
- + Smart search function
- + Direct connection to Bauerfeind BODYTRONIC measurement technology
- Patient file to store measurement and treatment data
- + Support via customer service, e-mail, contact form: experts.bauerfeind.com/contact

OUR COMPETENT BAUERFEIND PARTNERS

It is our aim also to ensure the quality of patient care, for the future as well. This can only be achieved by working with special suppliers who want to provide their customers with the best possible care as quickly as possible.

Our service concept has been approved precisely with the following in mind: efficient transactions, the best advice and tailor-made solutions for providing care.

Our sales staff can competently answer product-specific and business questions and are our contacts for joint marketing concepts. To make a rapid search or for convenient online ordering all the necessary information can be found at www.bauerfeind.de/find-bauerfeind.

SERVICE GUARANTEES **INTERNATIONAL** CUSTOMER SERVICE



100 % AVAILABILITY

By telephone within our service times (Mo. – Th. 7.00 am – 05.00 pm. Fr. 7.00 am – 04.00 pm CET). By fax, e-mail and electronic ordering systems 24 hours per day.



ELECTRONIC ORDERING SYSTEM

Make use of our offers for a fast, faultless and effort reduced data transfer: online shop, EDI service partner.



DESPATCH HANDLING

At your request we take care of the full processing of your shipment (incl. shipment guotations. additional documents) and we can bundle your orders flexible.



PROCESSING OF RETURNS

Within 5 days after receipt of goods.



ONE FACE TO THE CUSTOMER You have one fixed contact person for all of your desires and concerns.



PROACTIVE INFORMATION

We keep you posted – you receive from us all important information about your order and delivery status.



RETURN OF GOODS


Redemption of faultless goods within 30 days after delivery with 100 % refund or replacement (prior notification, filled form of return is required).

BAUERFEIND PRODUCT FEATURES A PLUS IN TERMS OF FUNCTION, QUALITY AND WEARING COMFORT

Bauerfeind has represented outstanding quality and excellent treatment for 90 years. Find out more about discerning expertise and high-quality product features that characterize our entire portfolio.

+ TRAIN KNIT FABRIC

The three-dimensional surface of the breathable material allows constant microcirculation and intermittent compression during movement. This helps to relieve pain and accelerates the reduction of edema. Muscular stabilization supports physiological movement processes. The Train knitted fabric fits perfectly and consists of different zones ensuring more effectiveness, comfort and temperature regulation.



+ COMFORT ZONE

For particularly sensitive areas of the body, such as crooks, bony projections and peripheral areas, we offer particularly stretchy material ensuring more wearing comfort, minimal pressure and easier doffing and donning.



+ GENTLE ON THE SKIN AND BREATHABLE

We really want to comply with the requests voiced by many users for particularly skin-friendly products. That is why we develop materials that are pleasant and breathable on the skin and fit without slipping.



+ EASY TO HANDLE

We want to make sure that our products are easy to put on and can be used without much effort. To achieve this, we integrate individual product solutions such as donning aids with a good grip, practical hand straps, special zippers, flexible Velcro straps or individual fastening systems.

+ TOOL-FREE FITTING

For the perfect adjustment to meet individual requirements, you can shape or configure many of our components, such as aluminum frames or hinges, without tools. Prompt service and uncomplicated use ensure a high degree of customer satisfaction and time savings.

+ MULTI-STAGE TREATMENT

The modular design of modern orthoses allows multi-stage treatment, depending on the required stabilization or desired mobilization by adding or removing functional elements.

+ PROVEN BY RESEARCH

Numerous clinical studies carried out by independent institutes show the medical effectiveness of the innovative materials, components and procedures of our products. Find out more about the latest study results from your Bauerfeind field service or in the certified CME training at our Academy.

+ DIGITAL MEASUREMENT

With the pioneering BODYTRONIC measurement system, Bauerfeind offers the option of individually adapting supports, medical compression stockings, foot orthoses and shoes, guaranteeing the perfect fit, excellent wearing comfort and maximum medical effectiveness – for highquality treatment.









THE NEW TRAIN SUPPORTS THE BENCHMARK FOR THE FUTURE

The therapeutic goals of medical supports are less pain and more joint stability.

But a support can only help when it fits properly and is comfortable to wear. It has to be convincing every day, after all:

THIS IS IMPORTANT TO PATIENTS WHEN BUYING A SUPPORT*

Bauerfeind has been acting according to these characteristics from the beginning.

The new Train supports are now even more

CONVINCING IN TERMS OF EFFECTIVENESS, WEARING COMFORT AND HANDLING.



to wear





easy and uncomplicated to put on and take off a sense of stability and confidence



REDUCING PAIN IN THE LONG TERM THE NEW TWO-COMPONENT PADS

- Pads provide relief and reduce pain
- Innovative massage nubs made of a second component now offer higher-intensity stimulation to ensure a long-lasting effect

MAKING DONNING AND DOFFING EVEN EASIER THE NEW PLASTIC STAYS

- Comfortable donning and doffing aid
- Durable, highly flexible synthetic material
- Adaptable flex zone ensuring the perfect fit

MORE INTENSE EFFECT THE NEW HIGH-LOW KNITTED FABRIC

+ STABILIZING, PAIN-REDUCING EFFECT

The highly structured surface creates a particularly stimulating massage and fits better than a smooth surface.

+ STIMULATES METABOLIC AND HEALING PROCESSES

The compressive knitted fabric creates a beneficial intermittent compression massage during movement.

+ EFFECTIVE THANKS TO THE PERFECT FIT

The knitted fabric is anatomically contoured and highly elastic to ensure a secure fit.

PLEASANT TO WEAR

The knitted fabric is thinner, the structure airier and compression more gentle in the sensitive areas of the body.

PROVEN BAUERFEIND QUALITY

COMFORTABLE AND DURABLE

Our elastic knitted material is surrounded by two soft yarns running in opposite directions to provide protection

SOFT EDGES

Wide edges with reduced pressure prevent pinching and guarantee unrestricted metabolism.



"THE FIT IS EXCEPTIONAL"

The three-dimensional, anatomical basic shape and high elasticity of the knitted fabric ensure that the new Train supports fit perfectly during movement, therefore being able to offer maximum effectiveness.

> Birgit Görnert, Head of Global Product Management


MEDICAL AID APP AND TREATMENT APP BAUERFEIND'S INTEGRATED CONCEPT

THE BAUERFEIND MEDICAL AID APP THE PRODUCT FINDER FOR PHYSICIANS AND MEDICAL RETAILERS

The Bauerfeind Medical Aid App is replacing the Physician App! Previously known as the Physician App, the digital product finder is now being launched with a new design.

This app will make choosing the right medical aid much easier. Thanks to its intuitive use, information and a current overview of our products are available everywhere – quickly and easily.

all Bauerfeind	11:00	100 % 🗰
-	BAUERFEIND"	
Q, Insti	atte oder Profilå) sycher	
BANDAGEN &	ORTHESEN	
1	1 6	-
FLB	Knie, Ober- schenkel & Hüfte	Isaule
1	scherikel & Hüfte	
(and	S) (1
Hand	Ellenbogen Sch	ulter
KOMPRESSIO	NSSTRÜMPFE	
EINLAGEN		
Haben Sie eine	Frage?	
§-Tipps		





Available free of charge for iOS and Android. For more information, please see: www.bauerfeind.de/himi-app

THE BAUERFEIND TREATMENT APP INDIVIDUAL TRAINING FOR YOUR CUSTOMERS



The body can only stay mobile and recover properly with training and movement. The Bauerfeind Treatment App supports patients to exercise with products from Bauerfeind.

Our trainer targets different problems and motivates users to do something to improve their health.

回認識個	
66.23	
10000	
间的外部	

Available free of charge for iOS and Android. For more information, please see: www.bauerfeind.de/app



THE TREATMENT APP IS A MEDICAL DEVICE.


Bauerfeind AG declares that the Treatment App is a Class I medical device in accordance with Appendix VII of the EU Medical Devices Directive 93/42/EEC. The app complies with the stipulations of the EU Medical Devices Directive 93/42/EEC.



DEVELOPED BY EXPERTS

The exercises have been compiled by physicians and physiotherapists to suit the relevant conditions.

WHAT CAN THE TREATMENT APP DO?





FOR TARGETED ANKLE TREATMENT

15

INDICATION GUIDE FOR FOOT AND ANKLE

PRODUCT		MalleoTrain®	MalleoTrain [®] S / S open heel	MalleoTrain° Plus	AchilloTrain°	AchilloTrain [®] Pro
			J.			
PRODUCT FEATURE	PRODUCT FEATURE		Stabilization of the ankle joint thanks to functional strap system	Increased stabilization and relief of the ankle thanks to functional strap system	Relief of the Achilles tendon thanks to heel wedge	Reduction of Achilles tendon problems thanks to friction pad
PAGE		18	19 + 20	21	22	23
INDICATIONS						
	mild	•	•	•		
Ligament insufficiency / instability	moderate					
	severe					
Irritation (post-traumatic/post-opera	tive)	•		•	•	•
Injuries to the capsular ligar	ment (acute)					
Early functional treatment		•		•		
Non-surgical/post-operative	e treatment	•		•		
Rehabilitation / protection		•	•	•		
Sinus tarsi syndrome / instability of the subtalar joi	nt					
Relapse prevention		•	•	٠		
Osteoarthritis/arthritis		•				
Tendomyopathy		•				
Achillodynia					•	•
Haglund's deformity					•	•
Hallux valgus (non-surgical	/ post-operative)					

AirLoc°	MalleoLoc [°] L	MalleoLoc° L3	MalleoLoc®	CaligaLoc®	ValguLoc°	ValguLoc [°] II
	I			Ot.	and a	
Stabilization of the upper ankle using adjustable air cushioning	Lateral stabilization of the ankle	Lateral stabilization of the ankle – removable in three stages	Increased stabilization of the upper and lower ankle	Stabilization of the upper and lower ankle after serious injuries	Correction of the big toe at rest	Correction of the big toe at rest and during movement
24	25	26	28	29	30	31
	•	•				
•	•	•	•			
			•	•		
	٠	•				
٠	•	•	•	•		
•	•	•	•			
•	•	•	•	•		
٠	•	•	•	•		•
				•		
•	•	•	•			
					•	•

MalleoTrain[®]

Relief and stabilization of the ankle









Lateral pad

Medial pad



- INDICATIONSChronic, post-operative and post-traumatic irritation
- (e.g. after sprains)Osteoarthritis and arthritis, joint effusions and swelling
- Tendomyopathies
- Ligament insufficiency
- Early functional/non-surgical treatment
- Relapse prevention

PROVEN IN A STUDY[△]

 80% of test subjects were pain-free after 14 days.

Pads anatomically contoured to the medial and lateral malleolus

BENEFITS AND MODES OF ACTION

- The regulated level of compression promotes the absorption of edemas, effusions and hematomas
- The pad redistributes support's compression away from the prominent malleoli to the surrounding soft tissues
- Neuromuscular joint stabilization by means of compressive knit material
- Intermittent compression is exerted during movement to promote circulation via the joint capsule and tendon compartments
- Optimized durability, even with intensive use, thanks to reinforced seams with an elastic band

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: right, left

Optional: Custom-made solution (beige, titan)

SCOPE OF DELIVERY

Including donning aid

SIZE	1	2	3	4	5	6	ORDER NO.
Ankle circumference (cm)	17–19	19-21	21-23	23-25	25-27	27-29	110 1110 20*0 +0°
* = Color (1) beige, (7) black, (+ = Design (6) right, (7) left ° = Size (1-6)	8) titan		points on	nade solution our order pa et for custor			

^a O'Hara J., Valle-Jones C. J., Walsh H, O'Hara H., Davey N. B., Hopkin-Richards H. and Butcher R. M. Controlled trial of an ankle support (MalleoTrain) in acute ankle injuries, Br J Sp Med 1992; 26(3)

experts.bauerfeind.com/malleotrain

MalleoTrain[®] S

Stabilization of the ankle joint thanks to functional strap system



Strap system wrapped in a figure of eight

Strap system

INDICATIONS

- Ligament insufficiency
- Supination prevention
- Mild sprains

BENEFITS AND MODES OF ACTION

- Stabilizes the ankle's movement planes and counteracts supination, thanks to its threelevel strap system (elastic, non-elastic, semi-elastic)
- Works just like a functional tape bandage, connecting the midfoot with the distal lower leg via a vertical figure of eight
- Optimized durability, even with intensive use, thanks to reinforced seams with an elastic band
- Neuromuscular joint stabilization using compressive knit material

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: black, titan Designs: right, left

SCOPE OF DELIVERY

Including donning aid

SIZE	1	2	3	4	5	6	ORDER NO.
Ankle circumference (cm)	17–19	19-21	21-23	23-25	25-27	27–29	110 1111 20*0 +0°
* = Color (7) black, (8) titan + = Design (6) right, (7) left							

° = Size (1 – 6)

FOOT & ANKLE

MalleoTrain[®] S open heel

Stabilization of the ankle joint and improved contact with the ground thanks to open heel









Strap system wrapped in a figure of eight

Strap system

Open heel

INDICATIONS

- Ligament insufficiency
- Supination prevention
- Mild sprains

BENEFITS AND MODES OF ACTION

- Stabilizes the ankle's movement planes and counteracts supination, thanks to its threelevel strap system (elastic, non-elastic, semi-elastic)
- Works just like a functional tape bandage, connecting the midfoot with the distal lower leg via a vertical figure of eight
- Neuromuscular joint stabilization by means of compressive knit material
- Optimized durability, even with intensive use, thanks to reinforced seams with an elastic band
- Open heel for improved contact with the ground

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, titan Designs: right, left

SIZE	1	2	3	4	5	ORDER NO.	
Ankle circumference (cm)	17–19	19-21	21-23	23-25	25-27	110 1115 00*0 +0°	~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~

* = Color (1) beige, (8) titan

+ = Design (6) right, (7) left ° = Size (1-5)

MalleoTrain[®] Plus

Increased stabilization and relief of the ankle thanks to functional strap system







8



Lateral pad

Medial pad

INDICATIONS

- Chronic, post-traumatic or post-operative irritation
- Post-operative rehabilitation
- Early functional treatment of capsular ligament injuries in the upper / lower ankle
- Ligament insufficiency / instability
- Supination prevention, particularly during increased activity

Strap system wrapped in a figure Strap system of eight

BENEFITS AND MODES OF ACTION

- Stabilizes the ankle's movement planes and counteracts supination, thanks to its threelevel strap system (elastic, non-elastic, semi-elastic)
- Works just like a functional tape bandage, connecting the midfoot with the distal lower leg via a vertical figure of eight
- Additional stability thanks to medial and lateral pads that are wrapped up tightly in a figure of eight pattern
- Neuromuscular joint stabilization by means of compressive knit material
- Optimized durability, even with intensive use, thanks to reinforced seams with an elastic band

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan **Designs:** right, left

SCOPE OF DELIVERY

Including donning aid

SIZE	1	2	3	4	5	6	ORDER NO.	
Ankle circumference (cm)	17–19	19-21	21-23	23-25	25-27	27-29	110 1112 0080 +0°	~

+ = Design (6) right, (7) left ° = Size (1-6) FOOT & ANKLE

AchilloTrain[®]

Relief of the Achilles tendon thanks to heel wedge







Pad running along

the Achilles tendon

Pad



Compensating cushion

INDICATIONS

Integrated heel wedge

- Achillodynia (tendinosis, paratendinitis, bursitis subachillea)
- Haglund's deformity
- Chronic, post-traumatic or postoperative irritation (e.g. in the case of Achilles tendon ruptures)

BENEFITS AND MODES OF ACTION

- The pad ensures a uniform application of pressure and exerts a massage effect during movement, improving the local metabolism and reducing edema
- The integrated tendon pad is viscoelastic and anatomically contoured
- The removable, integrated viscoelastic heel wedge (6 mm high) relieves the Achilles tendon
- Coordinated muscle control is assisted thanks to enhanced proprioception

APPLICATION INSTRUCTIONS

• When the heel wedge is inserted, the compensating cushion has to be worn on the other side.



- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: right, left

Accessories / spare parts: see page 106 and price list

SCOPE OF DELIVERY

Including donning aid

Including compensating cushion A separate heel cushion is also provided for the unaffected leg to offset the length difference.

SIZE	1	2	3	4	5	6	ORDER NO.
Ankle circumference (cm)	17 – 19	19-21	21-23	23-25	25-27	27-29	110 1101 30*0 +0°

* = Color (1) beige, (7) black, (8) titan

+ = Design (6) right, (7) left ° = Size (1-6)

experts.bauerfeind.com/achillotrain



AchilloTrain[®] Pro

Reduction of Achilles tendon problems thanks to friction pad







Friction pad

Friction pad running along the Achilles tendon

INDICATIONS

- Achillodynia (tendinosis, paratendinitis, bursitis subachillea)
- Haglund's deformity
- Chronic, post-traumatic or postoperative irritation (e.g. in the case of Achilles tendon ruptures)
- Inflammatory processes in the Achilles tendon

BENEFITS AND MODES OF ACTION

- The two wings of the pad extend to insertion of the calf muscle tendon and stimulate the upper third of the tendon and the nerve endings (proprioceptors) at the muscle-tendon junction
- The pad's silicone nubs point inward and provide a pain-relieving friction massage during movement. This improves proprioception and therefore the sequence of movement
- Vertical grooves aid lymph drainage, while the channel in the pad also guides the Achilles tendon in its anatomical path of movement

APPLICATION INSTRUCTIONS

• If needed, a heel wedge with a shock absorbing heel cushion is recommended

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: Identical design for right and left

SCOPE OF DELIVERY

Including donning aid

Ankle circumference (cm) 17-19 19-21 21-23 23-25 25-27 27-29 110 1105 10*0 00°	SIZE	1	2	3	4	5	6	ORDER NO.
Ankle circumierence (cm) 17–19 19–21 21–23 23–25 25–27 27–29 110 1105 10 0 00	Ankle circumference (cm)	17 – 19	19-21	21-23	23-25	25-27	27-29	110 1105 10*0 00°

* = Color (1) beige, (7) black, (8) titar ° = Size (1 - 6) **AirLoc**[®] Stabilization of the upper ankle using adjustable air cushioning







Adjustable air cushion

Elastic flexible shell

INDICATIONS

- Treatment of (acute) injuries to the capsular ligament complex
- Early functional / non-surgical treatment, e.g. after ankle sprains and ruptured ligaments
- Post-operative rehabilitation / protection, e.g. after suture repair / ligament reconstruction
- Chronic ligament insufficiency / instability
- Relapse prevention

BENEFITS AND MODES OF ACTION

- Stabilizes the lateral capsular ligament structure of the upper ankle, thanks to anatomically contoured plastic shells, thus preventing the foot from spraining (particularly when used in combination with shoes)
- The perfect fit and excellent adjustment to the respective degree of swelling, thanks to inflatable air cushions
- Allows ideal adaptation to the individual foot width and the degree of swelling using an elastic shell connection
- Flexible donning thanks to four adjustable Velcro straps

APPLICATION INSTRUCTIONS

• Individual adjustment of the air cushions using pump supplied

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Tool-free fitting
- + Treatment App

ORDERING OPTIONS

Colors: titan

Designs: right, left

Size: universal

Accessories / spare parts: see page 106 and price list

SCOPE OF DELIVERY

Including inflation bulb (pump)

PROVEN IN A STUDY^Δ

• 34% fewer complications in ankle treatment compared with taping

SIZE	ORDER NO.
Available in one universal size	 120 1303 1080 +00

+ = Design (6) right, (7) left

⁴ Lardenoye, S., Theunissen, Ed., Cleffken, B., Brink, P., R., G., de Bie, R., A., Poeze, M.:

The effect of faing versus semi-rigid bracing on patient outcome and satisfaction in ankle sprains: a prospective, randomized controlled trial BMC Musculoskeletal Disorders 2012, 13:8; http://www.biomedcentral.com/1471-2474/13/81

MalleoLoc[®] L

Lateral stabilization of the ankle







Microfiber cushioning material





Fits into everyday or business shoes

Plantar pad

INDICATIONS

- Non-surgical treatment of sprains and ruptured ligaments
- Acute and chronic capsular ligament instabilities
- Prevention of chronic ligament instability
- Post-operative rehabilitation / protection, e.g. after suture repair / ligament reconstruction
- Post-traumatic / post-operative irritation

BENEFITS AND MODES OF ACTION

Slim L shell

- Stabilizes the foot in the shoe and provides lateral protection against sprains, without restricting the foot's normal heel-to-toe movement
- The L shell sits in front of the talus bone and prevents talar shift
- The plantar pad (optional) that can be individually positioned stimulates the dorsiflexor muscles for additional active stabilization
- Excellent wearing comfort thanks to soft microfiber cushioning

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Tool-free fitting
- + Treatment App

ORDERING OPTIONS

Colors: titan

Designs: right, left

Size: universal

Accessories / spare parts: see page 106 and price list

PROVEN IN A STUDY^Δ

Counteracts spraining

SIZE ORDER NO. Available in one universal size 120 1326 1080 +00

+ = Design (6) right, (7) left

^a Gehring, D., Münch, M., Gollhofer A., Laborstudie: Nachweis von antisupinatorischen Effekten von Orthesen innerhalb des Verletzungsmechanismus; Orthopädieschuhtechnik, Ed. 10, p. 26–29, 2018

experts.bauerfeind.com/malleoloc-l

MalleoLoc[®] L3

Lateral stabilization of the ankle – removable in three stages





Treatment solution for three stability requirements

Zipper

Fits into everyday or business shoes

INDICATIONS

- Non-surgical treatment of sprains and ruptured ligaments
- Acute and chronic capsular ligament instabilities
- Prevention of chronic ligament instability
- Post-operative rehabilitation / protection, e.g. after suture repair / ligament reconstruction
- Post-traumatic / post-operative irritation

BENEFITS AND MODES OF ACTION

- The orthosis with removable elements supports patients throughout their treatment and offers a solution for three stability requirements according to the progress of healing
- The support's soft tissue compression helps to reduce swelling more quickly, it improves proprioception, and activates the muscles
- Provides lateral stabilization of the foot in the shoe and protects against sprains, without restricting the foot's normal heel-to-toe movement
- The L shell sits in front of the talus bone to counteract talar shift
- Zipper for easy donning of the support

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Multi-stage treatment
- + Treatment App

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 106 and price list

PROVEN IN A STUDY^Δ

Counteracts spraining

SIZE	1	2	3	4	5	6	ORDER NO.	
Ankle circumference (cm)	17 – 19	19-21	21-23	23-25	25-27	27–29	120 1327 1080 +0°	~~
+ = Design (6) right, (7) left ° = Size (1 - 6)								



^A Gehring, D., Münch, M., Gollhofer A., Laborstudie: Nachweis von antisupinatorischen Effekten von Orthesen innerhalb des Verletzungsmechanismus; Orthopädieschuhtechnik, Ed. 10, p. 26–29, 2018

experts.bauerfeind.com/malleoloc-l3

MalleoLoc[®]

Increased stabilization of the upper and lower ankle







Plantar guide

Anatomically shaped, bilateral stabilization element

INDICATIONS

- Early functional / non-surgical treatment for injuries to the capsular ligament complex and the bifurcate ligament, e.g. after sprains and ruptured ligaments
- Post-operative protection after suture repairs / ligament reconstruction
- Chronic ligament insufficiency
- Relapse prevention

BENEFITS AND MODES OF ACTION

- Stabilizes the foot and protects against sprains
- Allows both plantar flexion and dorsal extension, and prevents prevent talar shift, without restricting the "normal" gait
- Activates the dorsiflexor muscles using a plantar guide, thus supporting the stabilizing effect and supination protection offered by the orthosis

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Tool-free fitting
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 106 and price list

PROVEN IN A STUDY^Δ

Counteracts spraining

SIZE	1 (STANDARD)	2 (WIDE)	ORDER NO.	II
Heel width (cm)	< 6	> 6	120 1301 3080 +0°	
+ = Design (6) right, (7) ° = Size (1 - 2)	left			[]

° = Size (1 − 2)

^a Gehring et al.: Expecting ankle tilts and wearing an ankle brace influence joint control in an imitated ankle sprain mechanism during walking. Gait & Posture 2014; 39 (3), 894–898.

CaligaLoc[®]

Stabilization of the upper and lower ankle after serious injuries







Pronation wedge to relieve the lateral ligaments

Foot position without pronation wedge



- Non-surgical treatment of severe ankle injuries
- Post-operative protection after suture repairs / ligament reconstruction
- Chronic instability of the upper and / or lower ankle where surgery is contraindicated
- Post-traumatic sinus tarsi syndrome and decompensated instability of the subtalar joint, also as temporary stabilization aid

- BENEFITS AND MODES OF ACTION
- Prevents non-physiological stress situations in the upper and, to a certain extent, lower ankle

Foot position with pronation wedge

Medial view

• Raising of the hindfoot's outer edge and relief of the lateral ligaments thanks to integrated pronation wedge

APPLICATION INSTRUCTIONS

• Individual adjustment by cutting the cushion in the event of severe swelling thanks to cutouts in the area of the lateral malleolus

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Tool-free fitting
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan Designs: right, left

PROVEN IN A STUDY[△]

• Ready to work again 23.3 days earlier compared with plaster cast treatment

SIZE	1	2	3	ORDER NO.
Shoe sizes as a guide	34 - 38	39-43	44-48	120 1300 4080 +0°

+ = Design (6) right, (7) left ° = Size (1 - 3)

^a Zwipp et al, Die konservativ-funktionelle Behandlung des Knöchelbänderrisses hat sich bewährt: 5-Jahres-Ergebnisse, Med. Orth. Tech., Ausg. 114, S. 122–126, 1994 © Gentner Verlag Stuttgart

FOOT & ANKLE

ValguLoc[®]

Correction of the big toe at rest





Postural splint

INDICATION

• Hallux valgus (non-surgical and post-operative)

BENEFITS AND MODES OF ACTION

- Corrects the misalignment of the big toe and stretches shortened parts of the capsule and / or soft tissues
- Postoperative loss of correction (e.g. due to scarring) can be avoided, and the casing helps protect the big toe from certain mechanical forces
- Continuously adjustable correction pressure using Velcro fastening
- High degree of wearing comfort thanks to soft strap and cushioning materials

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Tool-free fitting
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 107 and price list

SIZE	1	2	3	ORDER NO.	
Shoe sizes as a guide	34 - 36	37-40	41-46	120 1310 2080 +0°	

+ = Design (6) right, (7) left ° = Size (1-3)

experts.bauerfeind.com/valguloc

ValguLoc[®] II

Correction of the big toe at rest and during movement







Multi-dimensional hinge to adjust the abduction and adduction angle

Complete mobility of the first metatarsophalangeal joint or immobilizing it with adjustable flexion and extension angles

INDICATIONS

- Post-operative protection / rehabilitation following hallux valgus surgery
- Non-surgical functional treatment of big toe misalignments with mobilization

BENEFITS AND MODES OF ACTION

- Individual configuration of the correction angle for the big toe thanks to hinge that can be adjusted three-dimensionally
- Can be adapted to the individual shape of the foot and locked in a range of positions if required for post-operative splinting / immobilization of the metatarsophalangeal joint of the toe
- Relieving pressure on the metatarsophalangeal joint of the toe using correction
- Fits into wide outdoor shoes thanks to low-profile design
- Enables natural heel-to-toe movement of the foot

THE ADVANTAGES OF BAUERFEIND

+ Easy to handle

ORDERING OPTIONS

Colors: titan

Designs: Identical design for right and left

Size: universal

Accessories / spare parts: see page 107 and price list

PROVEN IN A STUDY[△]

• 92% of patients rated the wearing comfort as good and very good.

SIZE	ORDER NO.
Available in one universal size	120 1316 0050 000

^a Nicht-interventionelle Studie 2012, Handhabung und Compliance der Hallux-valgus-Orthese ValguLoc II; internal data, Bauerfeind AG

experts.bauerfeind.com/valguloc-ii

FOOT & ANKLE



FOR TARGETED **KNEE TREATMENT**

Model wearing the GenuTrain® OA knee orthosis.

INDICATION GUIDE FOR KNEE, THIGH AND HIP

PRODUCT			GenuTrain*	GenuTrain [®] A3	GenuTrain [®] P3	GenuTrain [®] S	GenuTrain [®] S Pro	GenuPoint*	SofTec [®] Genu
					F				
PRODUCT FE	EATURE		Activation, relief and stabilization of the knee joint	Activation and stabilization for complex knee problems such as osteoarthritis of the knee	Centering and securing the patella using an adjustable corrective strap	Increased stabilization of the knee joint using joint splints on the sides	Increased stabilization of the knee joint using joint splints on the sides and limitation of range of motion	Targeted relief and guidance of the patellar tendon	Active and passive stabilization of the knee joint – ideal for long-term treatment
INDICATIONS	S	PAGE	36	38	39	40	41	42	43
Irritation (post-tr post-operative)	raumatic /		•	٠	٠			•	
Recurring joint e	effusion		•	•					
Patellar lateraliz	ation			•	•				
Patellofemoral p	ain syndron	ne			•			•	
	mild to m	oderate	•	٠	٠	•	•		
Instability	severe to	complex							•
	functional e.g. non-c								•
Meniscus injury			•				•		•
Collateral ligame	ent injury								٠
Ruptured anteric cruciate ligamen									•
Osteoarthritis / arthritis of the	mild to m	oderate	•	•		٠	•		
knee	moderate	to severe							
Condition after c	artilage trar	nsplant							
Patella fracture (non-surgical or	post-operat	ive)							
Genu recurvatun	n (with cros	s-strap)							
Corrective osteo	tomy								
Knee immobiliza	ation								
Muscle injury									
(Relapse) preven	ntion		•		•	•	٠		•
Osteoarthritis of	the hip								
Hip joint dislocat	tion								
Following endop revision procedu									
Muscular imbala lumbar, pelvis ar									
Post procedure muscular insufficiency, pain and stiffness									
Impingement of the hip joint (FAI)									
Post procedure i (after THR, after		gery)							
Feeling of instab									
Loosening of hip surgery is not po	prosthesis								

SecuTec [®] Genu	GenuTrain [®] OA	SofTec [®] 0A	SecuTec [®] 0A	MOS-Genu	GenuLoc®	CoxaTrain [®]	MyoTrain [®]	SofTec [®] Coxa
S	8		E					
Stabilization by limiting the ange of motion or complex snee injuries	Targeted relief and stabilization for more movement in cases of osteoarthritis of the knee	Relief for medial knee osteoarthritis using adjustable air bladders	Relief and stabilization for osteoarthritis of the knee with leg axis misalignment	Correction and stabilization after complex knee injuries or for corrective osteotomy	Immobilization of the knee joint – stretched or at a preset flexion angle	Stabilization and relief of pain in the hip joint	Adjustable compression for the treatment of thigh muscle injuries	Increased stabilization and guidance of the hip joint
44	46	48	49	50	51	52	54	55
								-
	•							
•				•				•
•	•		•	•				
•				•				
•			● in cases of leg axis misalignment	•	_			
	•							
	•	•	•	•				
•			•	•	•			
				•				
				•				
					•			
							•	
							•	•
						•		•
								•
								•
						•		
						•		
						•		
						•		
						٠		
						•		

GenuTrain®

Activation, relief and stabilization of the knee joint



37











Omega+ Pad

Train knit fabric

Integrated fitting and removal aid

Comfortable in the popliteal region Comfort sizes

INDICATIONS

- Tendinopathy, tendomyopathy, ligament insertion degeneration, meniscopathy, chronic, posttraumatic, or post-operative irritation
- Recurring joint effusion
- Osteoarthritis of the knee/arthritis (joint wear/joint inflammation)
- (Relapse) prevention
- Feeling of instability

PROVEN IN A STUDY[△]

patella

• 89% of subjects experience

• Measurably less strain on the

Proven reduction of pain when

greater knee stability

wearing the support

BENEFITS AND MODES OF ACTION

- Demonstrable relief and stabilization of the knee joint thanks to targeted activation of the muscles stabilizing the joint
- Massage of menisci through meniscus wings at the sides helps to reduce pain
- Supports the metabolism and the knee's selfhealing process thanks to the stimulating effect of the Hoffa spots on the infrapatellar (Hoffa's) fat pads
- Complex stimulation of mechanical sensors improves coordination, sensorimotor function and movement patterns
- Improvement of mobility by massaging tense muscles, tendons and fascia
- Accelerates the absorption of edemas and effusions and supports the healing process in case of irritation
- The breathable active zone and particularly soft comfort zones of the anatomically shaped Train knit fabric enable ideal wearing comfort with optimal medical effects

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: Identical design for right and left

Optional:

Siliconization: Standard sizes (black, titan), Comfort sizes (beige, black, titan)

Custom-made solution (beige, titan): siliconized edge, additional donning aid, velcro fastening

Accessories / spare parts:

Silicone dotted band (titan), see page 107 and price list

SIZE		-		2	-		-	-	7	ORDER NO.
Chandand	Circumference top (cm)	35-38	38-41	41 - 44	44-47	47-50	50-53	53-56	56-59	110 / 120 / 0*4 00%
Standard	Circumference bottom (cm)	25-28	28-31	31-34	34-37	37-40	40-43	43-46	46-49	110 4120 60 🛦 00*
Comfort	Circumference top (cm)				47 – 50	50-53	53-56	56-59	59-62	
Comfort	Circumference bottom (cm)				34-37	37-40				110 4120 60 🛦 08°

* = Color (1) beige, (7) black, (8) titan

Siliconization (0) without, (9) with (only titan and black, Comfort sizes also in beige) • = Size (0-7)

Custom-made solution: Please find the exact measuring points on our order pad for custom-made products or the order sheet for custom-made products on page 119.



experts.bauerfeind.com/genutrain





🙆 14 cm 🙂 12 cm

GenuTrain[®] A3

Activation and stabilization for complex knee problems such as osteoarthritis of the knee







Patella pad with friction zones

Integrated corrective strap for stabilizing the patella



- Osteoarthritis and arthritis of the knee, chronic, post-traumatic and post-operative irritation, ligament insertion degeneration and tendomyopathy
- Knee pain in cases of functional instability following muscular imbalance
- Recurring joint effusion
- (Relapse) prevention
- Patellar lateralization



Expansion zone in the lower leg area Siliconization

BENEFITS AND MODES OF ACTION

- · Pain reduction thanks to the targeted effect of the anatomically contoured pad with friction zones
- Activation of joint-stabilizing muscles and improved of coordination thanks to the stimulation of mechanosensors in the skin and the muscle-tendon junction
- Accelerated healing in cases of complex knee problems thanks to intermittent compression massage during movement; this promotes the reduction of edemas and the absorption of effusions
- Easier donning and doffing of the support based on the unique expansion zone in the lower leg area
- Maximum wearing comfort thanks to reduced pressure on the edges, the particularly soft area at the back of the knee and the moisture-regulating properties of the Train knitted fabric

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable

() 14 cm () 12 cm

- + Train knit fabric
- + Comfort zone
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, titan Designs: right, left

PROVEN IN A STUDY^Δ

• Pain-free walking: 2.4-fold increase of distance without pain

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference top (cm)	38-41	41 - 44	44-47	47-50	50-53	53-56	110 /105 10*0 . 00	-
Circumference bottom (cm)	28-31	31-34	34-37	37-40	40-43	43-46	110 4125 10*0 +0°	

* = Color (1) beige, (8) titan

+ = Design (6) right, (7) left ° = Size (1-6)

^a Reer et al.: The effect of a knee support in osteoarthritis; Orthopädie Technik, 8/2005

experts.bauerfeind.com/genutrain-a3

GenuTrain[®] P3

Centering and securing the patella with adjustable corrective strap









Patella pad

Corrective strap

Muscle-relaxing pad

INDICATIONS

- Patellar lateralization
- High-riding patella (patella alta)
- After lateral release surgery
- Knee joint irritation (patellofemoral pain syndrome, patellar chondropathy, chondromalacia, patella tip syndrome / jumper's knee, runner's knee)
- Feeling of instability

BENEFITS AND MODES OF ACTION

- Centering and protecting the patella in its physiological movement thanks to functionally shaped pad and individually adjustable corrective strap; It stops the patella from drifting laterally
- Friction points on the distal pad section massage the infrapatellar fat pad, thus relieving pressure and alleviating pain
- Relaxation of the iliotibial tract and reduced tension on the patella thanks to unique muscle-relaxing pad
- Improved proprioception thanks to pad wing with proximal extension as far as the vastus medialis muscle
- Accelerated absorption of edemas and effusions thanks to medically effective compression and intermittent compression massage during movement

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: right, left

Optional:

Siliconization (only titan), siliconized edge, custom-made solution (beige, titan)

SIZE	1	2	3	4	5	6	ORDER NO.	~
Circumference top (cm)	39-42	42-45	45-48	48-51	51 – 54	54-57		
Circumference bottom (cm)	28-31	31 – 34	34-37	37-40	40-43	43-46	110 4140 30 A +0 ²	(a) (b) (c) (c) (c) (c) (c) (c) (c) (c) (c) (c

* = Color (1) beige, (7) black, (8) titan

° = Size (1-6)

Custom-made solution: Please find the exact measuring points on our order pad for custom-made products or the order sheet for custom-made products on page 120.



GenuTrain[®] S

Increased stabilization of the knee joint using joint splints on the sides







Pad

Joint splints on the sides

INDICATIONS

- Joint instability (mild to moderate)
- Osteoarthritis of the knee
- Arthritis (e.g. rheumatoid arthritis)

BENEFITS AND MODES OF ACTION

- Joint splints on the sides relieve and guide the knee joint in combination with the circular strap system on the upper and lower leg
- Flat-fitting and individually adjustable tensioning straps with micro Velcro fastener
- High degree of wearing comfort thanks to breathable, moisture-absorbing Train knitted fabric with special area at the back of the knee and reduced pressure at the edges
- Supports neuromuscular joint stabilization by targeting proprioceptors
- Pain alleviation thanks to intermittent compression massage of the parapatellar soft tissue structures during movement owing to viscoelastic pad
- Accelerated absorption of edemas and effusions thanks to compression and massage effect of the elastic knitted fabric

APPLICATION INSTRUCTIONS

• The anatomically contoured joint splints are thermo-formable and can be adapted to the individual anatomy

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: right, left

Optional:

Siliconized edge, custom-made solution (titan)

Accessories / spare parts:

see page 107 and price list $% \left({{{\left[{{{\left[{{{\left[{{{\left[{{{c}}} \right]}}} \right]_{i}}} \right]}_{i}}}} \right]_{i}}} \right)$

SIZE	1	2	3	4	5	6	7	ORDER NO.
Circumference top (cm)	38-41	41 - 44						- 110 4130 40*0 +0°
Circumference bottom (cm)	28-31	31-34						- 110 4130 40'0 +0*
* = Color (1) beige, (7) black, (8) titan	Custom-	made solut	ion: Please	e find the e	xact meası	uring point	s on our	



+ = Design (6) right, (7) left ° = Size (1-7) Custom-made solution: Please find the exact measuring points on our order pad for custom-made products or the order sheet for custom-made products on page 121.

experts.bauerfeind.com/genutrain-s

GenuTrain[®] S Pro

Increased stabilization of the knee joint using joint splints on the sides and limitation of range of motion







Limitation of knee flexion

Pad

INDICATIONS

- Joint instability (mild to moderate)
- Osteoarthritis of the knee
- Arthritis (e.g. rheumatoid arthritis)
- Meniscus injuries

Limitation of knee extension

BENEFITS AND MODES OF ACTION

- Limiting the range of motion to support treatment using restrictable hinges
- Joint splints on the sides relieve and guide the knee joint in combination with the circular strap system on the upper and lower leg
- Supports neuromuscular joint stabilization by targeting proprioceptors
- Accelerated absorption of edemas and effusions thanks to compression and massage effect of the compression knit
- Flat-fitting and individually adjustable tensioning straps with micro Velcro fastener
- Pain alleviation thanks to intermittent compression massage of the parapatellar soft tissue structures during movement owing to viscoelastic pad

APPLICATION INSTRUCTIONS

- Flexion and extension of the knee joint can be limited independently of one another (flexion: 30°, 60°, 90°) (extension: 20°)
- The anatomically contoured joint splints are thermo-formable and can be adapted to the individual anatomy

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: beige, black, titan

Designs: right, left

Optional:

Siliconized edge, custom-made solution (titan)

SIZE	1	2	3	4	5	6	7	ORDER NO.
Circumference top (cm)				47-50				110 (125 00*0 . 00
Circumference bottom (cm)				37-40				110 4135 00*0 +0°

* = Color (1) beige, (7) black, (8) titan + = Design (6) right, (7) left ° = Size (1-7)

Custom-made solution: Please find the exact measuring points on our order pad for custom-made products or the order sheet for custom-made products on page 122.

experts.bauerfeind.com/genutrain-s-pro



A 14 cm

GenuPoint[®]

Targeted relief and guidance of the patellar tendon





Pad with frictional nubs

INDICATIONS

 Irritation in the patellar tendon area (patellofemoral pain syndrome, chondromalacia patellae, patella tip syndrome / jumper's knee)

BENEFITS AND MODES OF ACTION

- Guiding and relieving the patellar tendon using highly targeted compression
- The viscoelastic pad with special pressure points massages areas that are typically painful, alleviates symptoms and improves coordination
- Dynamic pressure adjustment and needsbased compression thanks to individually adjustable strap
- Maximum freedom of movement thanks to anatomic shape

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: black, titan

Designs: Identical design for right and left

PROVEN IN A STUDY^Δ

• Less pain and improved proprioception

SIZE	1	2	3	ORDER NO.
Circumference (cm)	28 – 30	30-35	35-44	120 4240 00*0 00°

* = Color (7) black, (8) titan

° = Size (1−3)

^a Astrid J. de Vries; Inge van den Akker-Scheek; Svenja L. Haak; Ron L. Diercks; Henk van der Worp; Johannes Zwerver Effect of a patellar strap on the joint position sense of the symptomatic knee in athletes with patellar tendinopathy Journal of Science and Medicine in Sport; 20, 11, S. 986-991, 2017

SofTec[®] Genu

Active and passive stabilization of the knee joint - ideal for long-term treatment





individual compromise axis of

• Rupture of the anterior and / or

posterior cruciate ligament

(ACL / PCL) for non-surgical

treatment, pre-operative and

long-term rehabilitation and in

cases of chronic insufficiencies

(traumatic, degenerative; e.g.

• Functional instabilities, e.g. lack

Collateral ligament injuries

of neuromuscular control,

• To restrict the range of motion

meniscal implantation)

PROVEN IN A STUDY^A

(e.g. after meniscus refixation or

 Mechanical and neuromuscular stabilization of the knee joint

"unhappy triad")

non-copers

Severe and / or complex instability

delayed post-operative treatment,

INDICATIONS

rotation

Stable, circular frictional connection

in accordance with the four-point

principle

Joint splints on the sides with restrictable joints

BENEFITS AND MODES OF ACTION

- Stabilization of the knee joint according to the four-point principle thanks to tensionstable and non-elastic knit sections, joint splints on the sides and Velcro straps
- The three-dimensional adjustment of the orthosis hinges to the individual compromise axis of rotation avoids harmful constraining forces on the knee joint
- Automatically determines the individual compromise axis of rotation
- Correct positioning of the SofTec Genu on the leg using patella cut-outs in the pad
- Stops the patella from drifting thanks to corrective strap integrated in the pad
- Long-lasting secure fit thanks to anatomically contoured Velcro tensioning straps

APPLICATION INSTRUCTIONS

• The hinges permit restriction of extension and flexion in 10° intervals



+ Easy to handle

Pad with integrated corrective strap

- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: black, titan Designs: right, left

Accessories / spare parts: see page 108 and price list

7 cm 🙃 15 cm

SIZE	1	2	3	4	5	6	ORDER NO.
Circumference top (cm)	38-42	41-45	44-49	48-53	51 – 57	54-60	
Circumference center (cm)	29-32	32-35	35-38	38-41	41-44	44-47	122 4600 10*0 +0°
Circumference bottom (cm)	28-32	31-35	34-38	36-40	39-43	42-46	
Hinge protection available as accessory							129 4601 1070 000

Hinge protection available as accessory

* = Color (7) black, (8) titan

+ = Design (6) right, (7) left


° = Size (1-6)

^a Schwameder, H.; Sell, S.; Strutzenberger, G.: Effect of two different functional braces on laxity and functional achievements in patients with anterior cruciateligament ruptures, submitted to: The American Journal of Sports Medicine


experts.bauerfeind.com/softec-genu

SecuTec[®] Genu

Stabilization by limiting range of motion for complex knee injuries













Easy to put on

Four-point principle in cases of ACL (blue arrows) or PCL injuries (yellow arrows)

Upper/lowerleg frames Hinge fastening of different sizes can be combined

BAUERFEIND

+ Easy to handle

+ Tool-free fitting

+ Digital measurement

ORDERING OPTIONS

Colors: blue, green

Designs: right, left

Versions: assembled,

partly assembled, construction kit

different sizes can be combined

(max. difference: two sizes);

Accessories / spare parts: see page 108 and price list

Delivered as assembly kit

Upper and lower leg sections of two

Calf strap adjustment

+ Gentle on the skin and breathable

THE ADVANTAGES OF

INDICATIONS

- Rupture of the anterior and / or posterior cruciate ligament (ACL/PCL)
- After ligament surgery / ligament reconstruction
- Severe and / or complex instability (traumatic, degenerative; e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion (e.g. after meniscus refixation or meniscal implantation)
- Patellar fracture (non-surgical and post-operative)

PROVEN IN A STUDY^Δ

• Especially secure fit during movement

- **BENEFITS AND MODES OF ACTION**
- Six straps in combination with a hard frame offer comprehensive stabilization following ACL, PCL or complex ligament injuries in accordance with the four-point principle
- Easily positioned at the front and fastened by the patient, reducing stress on the injured knee
- Tool-free limitation and immobilization of the joint; angle is clearly visible from the outside
- Anatomically correct and secure fit
- Extremely lightweight, slim and stable, thanks to specially processed light metal alloy frame
- High degree of rotational stability based on anatomically contoured frame design, including for the edges of the tibia

APPLICATION INSTRUCTIONS

- Flexion and extension can be limited independently of one another: Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90° Extension: 0°, 10°, 20°, 30°, 45°
- The joint can be immobilized in steps of 0°, 10°, 20°, 30°, 45°

SIZE	1	2	3	4	5	6	7
Circumference top (cm)	38-42	41 – 45	44-49	48-53	51 – 57	54-60	60-72
Circumference bottom (cm)	28-32	31-35	34-38	36-40	39-43	42-46	46-51

2

1 - 4

1

1-3



COMBINATIONS

° = Size (1-7)

Proximal frame size

Distal frame size	
**- Color (0E) blue (22) groop	

Color (05) blue, (23) green + = Design (6) right, (7) left

^ = Proximal frame size > = Distal frame size

3

1-5

	SecuTec Genu
4619 0**0 +0°	partly assembled
4647 0**0 +^>	Construction kit

45

4 5 2-6 3-7


6

4 – 7

	ORDER NO.	🙆 17 cm
	121 4614 0**0 +0°	SecuTec Genu
-		

GenuTrain[®] OA

Targeted relief and stabilization for more movement with osteoarthritis









Innovative relief system

Boa® Fit System

3D knitted mesh

Pressure-free bending zone

Large donning and removal aids

INDICATIONS

- Medial or lateral osteoarthritis of the knee
- Instability (mediolateral)
- Unicompartmental cartilage bone lesion
- Condition after cartilage surgery and cartilage growth for necessary unilateral postoperative relief
- Following unilateral meniscal lesions
- Ahlbäck's disease (osteonecrosis)
- Osteochondrosis dissecans
- Bone bruising (bone marrow edema of the medial femoral condyle / tibial plateau)

BENEFITS AND MODES OF ACTION

- Noticeable relief for pain on the side of the knee via the innovative relief and stabilization system its circumferential bands according to the three-point principle
- Simple adjustment of relief using the Boa Fit System dial – without further use of the Velcro fastenings
- More activity thanks to its light weight and secure position on the leg
- Flexible joint splint, figure-of-eight strap, and 3D knitted mesh address the sensorimotor function and promote natural movement patterns
- Comfortable to wear even over longer periods of time due to the breathable 3D knitted mesh and the pressure-free bending zone in the hollow of the knee
- Two easy-to-grip loops at the upper and lower edge of the orthosis make it easy to put on and take off

APPLICATION INSTRUCTIONS

- Tool-free reshaping of the joint splint for individual fitting to the patient's leg
- Easy handling of the Boa Fit System: – Open the dial to relax the bands and to put on or take off the orthosis
 - Close the dial and turn it to regulate pain relief

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Tool-free fitting
- + Treatment App

ORDERING OPTIONS

Colors: black-blue

Designs:

Right: for medial relief; or for lateral relief of the left leg Left: for medial relief; or for lateral relief of the right leg

Optional:

Limitation set in 3 sizes see page 108 and price list

$ \frac{\text{Circumference top (cm)}}{\text{Circumference bottom (cm)}} \frac{40-47}{28-34} \frac{47-52}{34-39} \frac{52-58}{39-43} \frac{57-65}{43-46} \frac{60-72}{46-50} $	SIZE	1	2	3	4	5	ORDER NO.
	Circumference top (cm)	40-47	47-52	52-58	57-65	60-72	101 /1/5 0070 .00
	Circumference bottom (cm)	28-34	34-39	39-43	43-46	46-50	121 4185 0070 +0*


+ = Design (6) right medial / left lateral, (7) left medial / right lateral ° = size (1-5)



SofTec[®] OA

Relief for medial knee osteoarthritis using adjustable air bladders







Relieves according to the three-point principle

Duo-centric hinge

INDICATION

• Medial knee osteoarthritis



- BENEFITS AND MODES OF ACTION
- The functional element with the tensionstable knitted fabric and strap system relieves the medial knee compartment according to the three-point principle
- The support effect accelerates the reduction of soft tissue irritation, edematous swellings and joint effusions
- Maximum freedom of movement for the muscles and, at the same time, a high degree of tension stability
- The valgus pressure can be individually adjusted at any time using integrated air bladders
- The correct fit of the orthosis is ensured by two pads at the upper and lower edges of the patella

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Digital measurement

ORDERING OPTIONS

Colors: black, titan

Designs: right, left

Accessories / spare parts: see page 108 and price list

SCOPE OF DELIVERY

Including inflation bulb (pump)

PROVEN IN A STUDY[△]

- 54 % less pain after wearing for six weeks
- 1 km average increase in patient mobility to 3.5 km

SIZE	1	2	3	4	5	6	7	ORDER NO.
Circumference top (cm)	41 - 45	44-49	48-53	51 – 57	54-60	57-63	60-66	
Circumference center (cm)	32-35	35-38	38-41	41 - 44	44-47	47 – 50	50-53	122 4680 00*0 +0°
Circumference bottom (cm)	31–35	34-38	36-40	39-43	42-46	45-49	48-52	-

* = Color (7) black, (8) titan

+ = Design (6) right, (7) left ° = Size (1-7)

🍐 von Herrwaarden R., Eine neuartige valgisierende Orthese zur Behandlung der medialen Gonarthrose..., MOT, (3) 2005

KNEE, THIGH & HIP

SecuTec[®] OA

Relief and stabilization for osteoarthritis of the knee with leg axis misalignment





relief, no narrowing of the knee



Adjustable calf shell


Different sizes can be combined to

provide the perfect fit

INDICATIONS

width

- Medial or lateral knee osteoarthritis / arthritis (moderate to severe)
- Condition after cartilage transplant
- After meniscus refixation
- Ruptured anterior / posterior cruciate ligament (ACL / PCL) with leg axis misalignment



BENEFITS AND MODES OF ACTION

- · Correction and relief according to the threepoint principle (varus and valgus setting)
- Stabilization according to the four-point principle (drawer effect; anterior/posterior)
- The perfect fit and pressure distribution thanks to adjustable calf shell and lower leg frame
- High level of wearing comfort thanks to lightweight aluminum frame
- · Easy to put on and take off thanks to quick-fit fasteners
- The anatomically contoured condylar gel cushions with insertable wedges ensure perfect adjustment

APPLICATION INSTRUCTIONS

- Individual fitting of the distal frame for varus or valgus adjustment (25° each) with 3 mm Allen key
- Flexion and extension can be limited independently of one another (flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°) (extension: 0°, 10°, 20°, 30°, 45°)

4

48-53

36-40

· Limiting hinge elements for restricting mobility or immobilization (to be ordered separately)

5

51-57

39-43

6

54-60

42-46

7

60-72

46-51



BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: blue

Designs: right, left

Versions:

Realignment orthosis, assembly kit

Upper and lower leg sections of two different sizes can be combined (max. difference: two sizes); Delivered as assembly kit

Accessories / spare parts: see page 109 and price list



							ORDER NO.	
Proximal frame size 1	2	3	4	5	6	7	121 4691 0050 +0°	SecuTec OA
Distal frame size 1–	3 1-4	1-5	2-6	3-7	4 – 7	5-7	121 4692 0050 +^>	Construction kit

ight, (7) left ° = Size (1 – 7)

Circumference top (cm)

Circumference bottom (cm)

SIZE

experts.bauerfeind.com/secutec-oa



3

44-49

34-38

2

41 - 45

31-35

1

38-42

28-32

> = Distal frame size

MOS-Genu

Correction and stabilization after complex knee injuries or corrective osteotomy









Long design

Adjustment in varus and valgus direction

Varus / valgus adjustment possible up to 12°

INDICATIONS

- Relief of the affected compartment in cases of varus / valgus knee osteoarthritis
- Rupture of the anterior and / or posterior cruciate ligament (ACL / PCL)
- Severe and / or complex instability (traumatic, degenerative; e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion (e.g. after meniscus refixation or meniscal implantation)
- Genu recurvatum (with cross-strap)

MOS-Genu long:

- Corrective osteotomies
- Complex ligament reconstructions
- Fractures located close to joints, implant surgery and protection
- Patellar fracture (non-surgical and post-operative)
- Condition after cartilage transplant

BENEFITS AND MODES OF ACTION

- The varus / valgus adjustment relieves the relevant knee compartment and / or the collateral ligament according to the three-point principle and corrects the leg axis alignment
- The knee orthosis is based on the well-known four-point principle and thus counteracts the anterior drawer effect
- The frame design stabilizes the knee joint against varus / valgus forces; the long version has a longer lever
- Hyperextension is prevented by using the popliteal strap (can be ordered as a spare part)
- The straps' large contact surface ensures excellent transmission of force and a secure fit
- The monocentric hinge with physiologically optimized pivot minimizes possible constraining forces exerted by the orthosis on the knee joint

APPLICATION INSTRUCTIONS

- Varus / valgus adjustment possible up to 12°
- The monocentric hinges can be restricted in 10° increments for flexion and extension.
- To determine the size of the MOS-Genu knee orthosis, the patient's knee should be flexed at 30°. Check the angle of flexion using the MOS Genu goniometer. Maximum M / L knee width is measured at the femoral condyles using the MOS-Genu orthometer (without applying pressure).

The orthosis size is read from the inspection window. If the measurement is on the border between two sizes, the degree of swelling or atrophy and the width of the other knee joint should be taken into consideration.

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: blue

Designs: short, long, right, left

Accessories / spare parts: see page 109 and price list

PROVEN IN A STUDY^Δ

• **30% less strain** on the medial side

SIZE	1	2	3	4	5	ORDER NO.
Max. Knee width (mm)	85 – 95	95–105	105–115	115–125	125–135	121 4660 405# +0°



= Design (1) short, (2) long + = Design (6) right, (7) left

° = size (1-5)

experts.bauerfeind.com/mos-genu

Kutzner et al., The effect of valgus braces on medial compartment load of the knee joint – in vivo load measurements in three subjects, Journal of Biomechanics 44, 1354–1360, 2011

GenuLoc®

Immobilization of the knee joint – stretched or at a preset flexion angle











Plastic stays at the sides

Pad behind the knee

Contoured aluminum splints

INDICATIONS

- Pre-operative, post-operative, post-traumatic immobilization
- Non-surgical treatment of patella fractures

BENEFITS AND MODES OF ACTION

- The knee's flexion angle can be adjusted using dorsal aluminum splints
- Reduced edema development thanks to removable pad behind the knee
- Stabilizing plastic stays at the sides automatically adjust to the set knee flexion angle, guaranteeing an excellent fit

APPLICATION INSTRUCTIONS

- The aluminum splints are contoured to ensure a knee flexion angle of 20°
- Universal size, one design for the right and left knee

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable

ORDERING OPTIONS

Colors: black **Size:** universal

 SIZE
 ORDER NO.

 Available in one universal size
 120 4320 0070 000

CoxaTrain[®]

Stabilization and relief of pain in the hip joint













KNEE, THIGH & HIP

53

Stabilization

Gluteal pad – massages trigger points

hip joint

Tool-free adjustment

Height-adjustable joint splint

INDICATIONS

- Moderate to severe osteoarthritis (Kellgren-Lawrence score 2-4)
- Muscular imbalance in the lumbar, pelvis and hip area
- Post procedure muscular insufficiency, pain and stiffness
- Impingement of the hip joint (FAI)
- Post procedure instability (after THR, after revision surgery)
- Feeling of instability of the hip
- Loosening of hip prosthesis if surgery is not possible

BENEFITS AND MODES OF ACTION

- Pain-relieving counteraction of muscular imbalances in the lumbar, pelvis and hip area during movement with a positive effect on gait:
 - Stabilizes the pelvis and relieves the sacroiliac joints (SI joints)
 - Dorsal SI joint friction pads massage the local ligament structures and adjacent muscles
 - Trochanter pad moves during walking, providing friction massage
- Gluteal pad massages trigger points
- Stabilizing guidance provided by the lateral hinge
- High level of patient compliance
- Easy to put on without much effort thanks to looped tensioning straps and practical finger pockets/loops

APPLICATION INSTRUCTIONS

Customizable fitting:

- Tool-free adjustment of the hinge to suit the patient's anatomy
- Height-adjustable thigh support above the joint splint
- Individually adjustable tensioning straps on the pelvic support for stabilization

• Three degrees of restriction of hip movement (flexion and extension) thanks to optional joint limiting pin (60°, 75°, 90°)

ADVANTAGES OF BAUERFEIND

- + Easy to adjust
- + Gentle on the skin and breathable
- + Tool-free fitting
- + Latex free

ORDERING OPTIONS

Colors: titan (grey / blue)

Optional:

Short hinge < 170 cm body height Long hinge > 170 cm body height

Note:

• In specific cases, it is possible that a "long" hinge will fit better than a measured size of "short" or the other way round.

Accessories / spare parts:

see page 109 and price list

PELVIC SUPPORT SIZE	1	2	3	4	5	6	ORDER NO.
Circumference of pelvis (cm)	85-95	<u>-</u> 95–105	105–115	4 115–125	125-135	135-145	
° = size (1-6)							
THIGH SUPPORT							
SIZE	1	2	3				ORDER NO.
Circumference of thigh (center) (cm)	43-53	53-63	63-73				122 2582 1080 +0°
+ = Design (6) right, (7) left ° = size (1 – 3)							
HINGE							
SIZE	SHORT	LONG					ORDER NO.
Height (cm)	<170 cm	> 170 cm					122 2584 100# 000

= Design (1) short, (2) long

experts.bauerfeind.com/coxatrain

MyoTrain[®]

Adjustable compression for the treatment of thigh muscle injuries





Optional: functional pad and small acute pad

INDICATIONS

- For compression in the case of muscle injuries (muscle hardening, pulled muscles, torn muscle fibers, torn muscles, muscle bruising)
- Relapse prevention

BENEFITS AND MODES OF ACTION

- Works according to the principles of muscle suspension bands with compression in the injured area
- Regulation of compression using Velcro strap system and two freely placeable pads
- Targeted effect thanks to acute pad straight after injury to inhibit bleeding and limit hematoma
- Extensive compression thanks to large functional pad, minimizing swelling and pain
- Maximum wearing comfort and a secure fit thanks to anatomical knitted fabric as well as siliconized and pressure-reduced edges
- Regulates muscle tone and prevents recurring injuries

THE ADVANTAGES OF BAUERFEIND

+ Easy to handle

ORDERING OPTIONS

Colors: titan

Designs: Identical design for right and left

SIZE	1	2	3	4	5	6	7	ORDER NO.	~
Circumference top (cm)	44 - 48	48-52	52-56	56-60	60-64	64-68	68-72	- 110 3250 0080 00°	S @
Circumference bottom (cm)	32 – 35	35-38	38-41	41-44	44-47	47 – 50			30°

• = Size (1-7)

SofTec[®] Coxa

Increased stabilization and guidance of the hip joint







Monocentric hinge for abduction / adduction and flexion/extension

INDICATIONS

• Dislocation

surgery

procedures

Instability





Treatment for both sides

- Adduction and abduction can be adjusted continuously between -6° and +6°.
- The hinge is available in offset and straight versions and can be used for either the right or left hip
- The orthosis can be adapted directly at the hospital bed, even if the patient is lying down.

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Optional:

ORDER NO.

122 252# 0080 00°

Offset hinge / straight hinge

Accessories / spare parts: see page 110 and price list



Total replacement surgery

• Total replacement revision

(Girdlestone arthroplasty)

• Hip spacer for two-stage revision

• Femoral head resection

· Osteoarthritis of the hip

(moderate to severe)

Three-part plastic pelvic harness

BENEFITS AND MODES OF ACTION

- The orthosis holds the ball head securely in the socket using a combination of
 - secure positioning of the pelvic frame on the iliac crests
 - individual adjustment of the thigh harness on the sides
 - Adapting the orthosis hinge to suit the indication
- The multi-dimensional, monocentric hinge allows for treatment in the adduction or abduction position as required and limits extension and flexion
- Restricted, guided movement safeguards the outcome of surgery
- The anatomically shaped pelvic frame that is integrated into the knitted fabric guarantees that the orthosis can be positioned easily and correctly
- The knitted fabric sections and the tensioning straps on the pelvic frame and thigh harness can be trimmed to fit
- Perfect guidance of the hip joint and anatomically correct fit thanks to monocentric orthosis hinge

APPLICATION INSTRUCTIONS

3

115-135

2

95-115

- Can be adjusted individually and parallel to the hip's pivot point
- The extension and flexion of the hinge can be restricted in 10° increments from -10° to +90°

= Design (4) Offset hinge (for women), (5) straight hinge (for men) = Size (1 - 3)

1

75-95

Circumference trunk (cm)

SIZE







FOR TARGETED BACK TREATMENT

INDICATION GUIDE FOR SPINE AND NECK

PRODUCT		LumboTrain°/ Lady	SacroLoc®	LordoLoc*	LumboLoc*	LumboLoc° Forte	Spinova° Stabi Classic	Spinova° Support Plus	Spinova® Support Plus Classic	Spinova® Unload Classic
								5		
PRODUCT FEATUR	E	Activation, relief and stabilization of the lumbar spine	Stabilization of the pelvis and targeted relief of the sacroiliac joints	Stabilization and relief of the lumbar spine using tensioning straps	Stabilization and relief of the lumbar spine	Increased stabilization and relief of the lumbar spine using adjustable strap system	Moderate stabilization of the lumbar spine using corset stays	Stabilization of the lumbar spine using anatomically contoured back pad	Moderate support of natural lumbar lordosis using an anatomically shaped back pad	Moderate lordosis reduc- tion and relief of the lumbar spine using removable bridging frame
FUNCTIONAL ELEM	IENT		db		\ /	O or			removable	removable
PAGE		60	62	64	65	66	70	71	72	73
INDICATIONS										
Low back pain (acute and o	chronic)	•		•	•					
(Pseudo) radicular lumbar syndrome / lumbosciatica with muscular deficiencies	moderate					•	•	•	•	•
in cases of intervertebral disk protrusions	severe									
Degeneration /	mild	•								
muscular insufficiency of the spine	moderate			•	•	•			•	
	severe							•		
Condition following herniation (non-surgical)		•					•	•	•	
Condition following discert herniation (post-operative)		•				•	•	•	•	•
Facet syndrome / osteoarthritis	moderate 				•			•	•	•
Lumbar deformities						•		•		
Ligamento-tendomyopathi	es					٠				
	Grade I					•	•			•
Spondylolysis / spondylolisthesis	Grade II							•	•	•
	Grade III									
Spinal fusion / kyphoplasty (post-operative)	,							•	•	
Spinal canal stenosis								•	•	•
Foraminal stenosis										
Osteoporosis	al state							•	•	
SI joint syndrome, osteoar instability, blockage	thritis,		•							
Pelvic girdle instability, e.g rupture and loosening	ı. symphysis		•							
Myalgia, tendinopathy in th	ne pelvic region		•							
Structural disturbance foll fusion (L5/S1)	owing spinal		•							
Non-surgical treatment following pelvic fractures			•							
Vertebral Anterior ed fractures Posterior e	-							•	•	•
Osteochondrosis										
Tumors (metastases)										
Vertebral fractures (lumba	r spine)									
Painful hunched back										
Pain in the cervical spine	moderate									
	severe									
Whiplash / blockage										

Spinova® Unload Plus	Spinova° Immo Plus	Spinova [®] Immo Plus Classic	Spinova° Immo	SecuTec* Dorso	SecuTec° Lumbo	SofTec° Lumbo	DorsoTrain*	Spinova° Osteo	SofTec* Dorso	CerviLoc®	CerviLoc [®] S
	8		4				R				2
Lordosis reduction and increased relief of the lumbar spine using removable bridging frame	Immobilization of the lumbar spine thanks to mobilization function using removable shell	Immobilization of the lumbar spine thanks to mobilization function and removable shell	Immobilization of the lumbar spine	Lordosis reduction and relief of the lumbar spine using a bridging frame and corset stays	Lordosis reduction and relief of the lumbar spine using bridging frame	Immobilization of the lumbar spine thanks to mobilization function using removable shell	Active straightening and support of the spine	Stabilization and active straightening of the spine, with free abdomen and chest area	Increased stabilization of the thoracic and lumbar spine	Stabilization the cervical s	
J	removable	removable		removable	removable	removable	ſ		Ĭ		
74	75	76	77	78	79	80	81	82	84	85	85
•	•	•		•	•	•					
							•				
	•	•			•	•		•	•		
	•	•									
•	•	•		•	•	•			•		
•	•	•		•	•	•					
•	•			•	٠						
•	•			•	•						
•	•	•		•	•	•					
	•	•			•	•			•		
•	•	•		•	٠	•					
•	•	•		• 1			•	٠	•		
	•	•	•	• 1		•					
•	•	•	•	•	٠	•					
	•	•				•					
				•			• ²	• ²	• 3		
										•	
										•	•
1 without frame	² osteoporotic	³ single up to TI	-18							•	•

SPINE & NECK

LumboTrain[®] / Lady

Activation, relief and stabilization of the lumbar spine







Pad can be positioned individually





Practical hand straps make the support easy to put on and take off

Perforated closure and large interstices in the knit fabric for high breathability

LumboTrain Lady (fitted cut)

INDICATIONS

Pad with frictional nubs

- General pain in the lumbar spine region (acute and chronic)
- Degenerative conditions of the lumbar spine, e.g. osteochondrosis, osteoarthritis, lumbar deformities, spondylolysis without vertebral displacement
- Condition following herniation (non-surgical / post-operative)
- Condition following discectomy

BENEFITS AND MODES OF ACTION

- Targeted pain relief during movement is enabled by intermittent compression massage from the knit fabric and viscoelastic pad
- Easy to put on with minimal effort using practical hand straps
- Stabilizes and relieves the lumbar region thanks to the straightening effect of circular compression
- LumboTrain Lady with fitted shape and high cut leg

APPLICATION INSTRUCTIONS

 Pad can be positioned individually as required

BAUERFEIND + Easy to handle

THE ADVANTAGES OF

- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

Designs: straight / fitted cut

Accessories / spare parts: see page 110 and price list

SPINE & NECK

PROVEN IN A STUDY^Δ

• 21 % higher back muscle activity

• **69 % less pain** compared with the control group

SIZE	1	2	3	4	5	6	7	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-100	100-110	110-120	120-130	130-145	110 841# 5080 00°	- (

= Design (0) straight cut, (1) fitted cut
° = Size (1 - 7)

Anders, C., Hübner, A.

Plues 9, c. (nounce, A. Influence of elastic lumbar support belts on trunk muscle function in patients with nonspecific acute lumbar back pain PLoS ONE 14(1): e0211042. https://doi.org/10.1371/journal.pone.0211042; 2019

SacroLoc®

Stabilization of the pelvis and targeted relief of the sacroiliac joints













63

Two massage pads

Stabilization

Elastic mesh material with patented strap system

Individually adjustable tensioning straps

INDICATIONS

- Sacroiliac joint disorders (SI joint syndrome, osteoarthritis of the SI joint, SI joint instability)
- Pelvic girdle instability, e.g. symphysis rupture and loosening
- Sacroiliac joint blockage (acute and prophylactic)
- Structural disturbances following spinal fusion (in the transition from L5 to S1)
- Myalgia and tendinopathy in the pelvic region
- Non-surgical treatment following pelvic fractures

BENEFITS AND MODES OF ACTION

- Provides immediate relief to the pelvis and sacroiliac joints thanks to circular compression using individually adjustable tensioning straps
- Targeted pain reduction thanks to two friction pads on the sacroiliac joints which can be positioned precisely
- Exceptional level of wearing comfort, even when sitting – thanks to breathable mesh materials and narrow, low-profile fastening
- Deflection loops at the sides for a high level of pressure without exerting much energy
- Practical, wide hand straps making it convenient and easy to put on with minimal effort
- High level of patient compliance as the orthosis can be worn inconspicuously under clothing
- The rigid tensioning strap system ensures circular stabilization - nutation of the sacral bone is restricted, ligaments are relieved and pain is reduced

APPLICATION INSTRUCTIONS

- The suitability and use of the orthosis during pregnancy must be determined by a physician on a case-by-case basis
- Pads can be positioned individually as required

+ Digital measurement

+ Gentle on the skin and breathable

THE ADVANTAGES OF

+ Treatment App

BAUERFEIND

+ Easy to handle

ORDERING OPTIONS

Accessories / spare parts:

PROVEN IN A STUDY^A

• 32% pain reduction

SIZE	1	2	3	4	5	6	7	ORDER NO.
Circumference pelvis (cm)	70 - 80	80-90	90-100	100-110	110-120	120-130	130-140	120 8416 2080 00°
° = Size (1 – 7)								



^a Sichting et al. Pelvic Belt Effects on Sacroiliac Joint Ligaments: A Computational Approach to Understand Therapeutic Effects of Pelvic Belts Pain Physician 2014; 17:43-51 · ISSN 1533-3159. Hammer et al. Pelvic Belt Effects on Health Outcomes and Functional Parameters of Patients with Sacroiliac Joint Pain PLOS ONE | DOI:10.1371/journal.pone.0136375 August 25, 2015

Colors: titan

see page 110 and price list

LordoLoc[®]

Stabilization and relief of the lumbar spine using tensioning straps







Integrated dorsal corset stays

Elastic tensioning straps for Stabilization and relief individually adjustable compression

INDICATIONS

- General pain in the lumbar spine region (acute and chronic)
- Moderate degeneration / muscular insufficiency



- Integrated corset stays straighten the lumbar spine
- Increased stabilization thanks to individually adjustable tensioning straps
- Unobtrusive under clothing thanks to lightweight, low-profile design
- The compressive, elastic material relieves and activates the abdominal and back muscles

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	70-80	80-90	90-100	100-110	110-120	120-130	121 8500 0080 00°	



LumboLoc[®]

Stabilization and relief of the lumbar spine







Integrated dorsal corset stays straighten the lumbar spine.

Stabilization and relief

INDICATIONS

- General pain in the lumbar spine region (acute and chronic)
- Moderate degeneration / muscular insufficiency
- Facet syndrome

BENEFITS AND MODES OF ACTION

- Integrated corset stays straighten the lumbar spine
- Stabilizes and relieves the lumbar spine thanks to the compressive Train knit
- Easy to put on using practical hand straps
- The perfect fit based on elastic material and adaptable corset stays

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Train knit fabric
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

SPINE & NECK

PROVEN IN A STUDY^Δ

• 53 fewer painful days per year



^a Roelofs et al. Lumbar Supports to Prevent Recurrent Low Back Pain among Home Care Workers Annals of Internal Medicine 2007, 147, 685-692

experts.bauerfeind.com/lumboloc

LumboLoc[®] Forte

Increased stabilization and relief of the lumbar spine using adjustable strap system







D



Sacral pad

Lumbar pad

Corset stays

Functional strap

Stabilization and relief

INDICATIONS

- Moderate lumbar sciatica with slight muscular deficiencies in cases of intervertebral disk protrusions / herniated disks
- Moderate (pseudo) radicular lumbar syndrome
- Spondylolysis / spondylolisthesis, grade 1
- Moderate lumbar deformity
- Condition following herniation, e.g. post-operative
- Ligamento-tendomyopathy
- Moderate degeneration / muscular insufficiency

BENEFITS AND MODES OF ACTION

- Increased stabilization thanks to individually adjustable tensioning straps
- Integrated corset stays straighten the lumbar spine
- Height-adjustable functional strap for indication-specific pressure distribution
- Additional stabilization using pads for the lumbar area or sacral bone*

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

Optional: Indicate optional lumbar or sacral pad in prescription; must be ordered separately

Accessories / spare parts:

see page 110 and price list

* Please include the type of pad on the prescription.

LumboLoc Forte

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	70-80	80-90	90–100	100 – 110	110-120	120–130	121 8403 0080 00°	\sim

 $^{\circ}$ = Size (1-6)

LumboLoc Forte pads

SIZE	ORDER NO.
Available in one universal size	121 840# 0050 001
# = Design (4) sacral pad, (5) lumbar pad	

experts.bauerfeind.com/lumboloc-forte

Spinova[®]-System

The modular Spinova system for targeted treatment in cases of lumbar spine indications

The Spinova system with six different back orthoses is an innovative product concept by Bauerfeind, designed to enhance the range of treatment options on offer for various indications relating to the lumbar spine.

The Spinova back orthoses are based on the same concept: the comfortable Spinova mesh material is the general functional element. Thanks to integrated circular corset stays, it acts like a Lindemann corset. Circular force is exerted where relevant to the treatment, using the adjustable Spinova tensioning strap system with back setting device. The Classic Line's tensioning strap system exerts moderate force. The system is supplemented by indication-specific Spinova functional elements, such as a back or abdominal pad, bridging frame or shell.

Spinova lumbar orthoses are used to correct posture or to provide relief to the lumbar spine. They are an effective method for preventing non-physiological strain on the spine and the intervertebral disks.

They can be adapted in many ways and provide the patient with unique wearing comfort and ensure easy handling.

CLASSIC TENSIONING STRAP SYSTEM

for mild to moderate indications

- Moderate level of pressure
- Individually adjustable compression
- Convenient and easy to put on with minimal effort

Product variants:

- Spinova[®] Stabi Classic
- Spinova[®] Support Plus Classic
- Spinova[®] Unload Classic
- Spinova[®] Immo Plus Classic

ADJUSTABLE TENSIONING STRAP SYSTEM

for moderate to severe indications

- Perfect stability thanks to targeted pressure distribution
- Indication-specific adjustment
- Direction and intensity can be adapted
- Can be safely reproduced by patient

Product variants:

- Spinova[®] Support Plus
- Spinova[®] Unload Plus
- Spinova[®] Immo Plus
- Spinova[®] Immo






Spinova[®] system



Adjustable tensioning strap system, front view



Adjustable tensioning strap system, rear view



reduction

enhancement

GENERAL FUNCTIONAL ELEMENTS: SPINOVA SYSTEM



Integrated circular corset stays corresponding to Lindemann corset/corset based on a grid of stays



Elastic Spinova mesh material, anatomically contoured, highly breathable



Anatomically shaped abdominal fastening with micro Velcro and finger pockets for easy donning and doffing

FUNCTIONAL ELEMENT CLASSIC **TENSIONING STRAP SYSTEM**



Classic Velcro straps that can be individually tightened provide the lumbar spine with moderate stabilization

FUNCTIONAL ELEMENTS ADJUST-ABLE TENSIONING STRAP SYSTEM



Adaptable back setting device stabilizes and corrects the area from the sacral bone to the lower thoracic spine



The back setting device positioned parallel to the spine bridges spinous processes; no pressure on the wound

Individually heightadjustable tensioning straps ensure indicationspecific pressure distribution



Strong, defined tension thanks to adjustable torque level in the strap closure for stabilization



Patented Spinova tensioning strap system ensures tensile force from three directions



Consistent therapeutic circumferential pressure thanks to adjustable strap tension and buckle fastening

INDICATION-SPECIFIC MULTI-STAGE TREATMENT CONCEPT

The modular design of the Spinova lumbar orthoses allows individual treatment by adding or removing indication-specific functional elements in accordance with a two-stage treatment process.

INDICATION-SPECIFIC FUNCTIONAL ELEMENTS CAN BE ADDED AND REMOVED

Bridging frame

Pads



- Anatomically self-adjusting, pivoting cross stays
- Easily malleable, heightadjustable aluminum stays



- Back pad ensures extensive bridging and support for lumbar lordosis
- Abdominal pad offers increased compression of the abdomen

Shell system



- Individually adjustable thanks to shell sections that can be moved in relation to each other
- Shell can be removed without tools by loosening the elastic bands and knurled screw
- Can be used immediately thanks to preset shell circumference
- Less adjustment required thanks to two shell sizes (short: 33 cm, standard: 38 cm)

TWO-STAGE TREATMENT CONCEPT

1st stage: stabilization phase

For the effective post-operative or initially conservative correction of segments, the orthosis is used together with the indication-specific functional elements (shell / bridging frame / pads).

2nd stage: mobilization phase

When the patient becomes more mobile, the indicationspecific functional elements are removed. The corset stays that remain in the orthosis provide external support.

Elements can be added to the orthosis at any time if required, depending on the patient's stabilization needs.



69

SPINE 8 NECK

FIND OUT MORE AT https://player.vimeo<u>.com/</u>

video/158159586

Spinova[®] Stabi Classic

Moderate stabilization of the lumbar spine using corset stays





Integrated circular corset stays, can be adjusted

Finger pockets on the abdominal fastening

INDICATIONS

- Moderate (pseudo) radicular lumbar syndrome
- Moderate lumbosciatica with slight deficiencies in the case of intervertebral disk protrusions / herniated disks
- Condition following herniation (non-surgical / post-operative)
- Condition following discectomy
- Spondylolisthesis, grade I

BENEFITS AND MODES OF ACTION

- Moderately stabilizes and straightens the lumbar spine thanks to integrated circular corset stays
- The supportive effect extends as far as the lower thoracic spine
- Easy to put on using large finger pockets on the abdominal fastening
- The perfect fit thanks to adaptable aluminum stays
- The length of the tensioning straps can be individually adjusted

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement
- + Treatment App

ORDERING OPTIONS

Colors: titan

Accessories / spare parts: see page 112 and price list

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-105	105-115	115–125	125–140	122 8561 0000 00°	
° = Size (1-6)								

Spinova[®] Support Plus

Stabilization of the lumbar spine using anatomically contoured back pad





Integrated circular corset stays,



strap system ensuring regulated pressure distribution



Lordosis enhancement

INDICATIONS

can be adjusted

- Vertebral fractures (lumbar spine) with vertebral body damage relieving the anterior edges of vertebrae
- Spondylolysis / spondylolisthesis, grade II
- Moderate lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Lumbar spinal canal stenosis, also with paresis (non-surgical, post-operative)
- Moderate (pseudo) radicular lumbar syndrome
- Severe degeneration / muscular insufficiency
- Condition following herniation (post-operative / non-surgical)
- Condition following discectomy
- Condition following spinal fusion / kyphoplasty (one or more levels)
- Osteoporosis

PROVEN IN A STUDY^Δ

• Good patient compliance and great sense of stability

В	ENEFITS AND MODES OF ACTION
•	Integrated circular corset stays ensure

- Int stabilization and straightening
- Individually adaptable tensioning strap system for regulated pressure distribution
- Greater straightening thanks to increased abdominal pressure exerted by the abdominal pad
- Supports natural lumbar lordosis using a large, anatomically shaped back pad

THE ADVANTAGES OF BAUERFEIND

- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, standard

- short: body height < 170 cm
- standard: body height > 170 cm

Accessories / spare parts: see page 112 and price list

SIZE	1	2	3	4	5	6	7	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-100	100-110	110-120	120-130	130-140	122 8542 100+ 00°	<
+ = Design (0) standard, (1) sho	ort								

• = Size (1-7)

Rommelspacher, Y., Ziob, J., Wirksamkeit der Rückenorthese Spinova Support Plus bei der postoperativen Versorgung von segmentalen Spondylodesen degenerativer LWS-Erkrankungen, Department for Orthopedics and Trauma Surgery, Bonn University Hospital, submitted for publishing

experts.bauerfeind.com/spinova-support-plus



Abdominal pad

Spinova[®] Support Plus Classic

Moderate support of natural lumbar lordosis using an anatomically shaped back pad





Removable back pad

INDICATIONS

- Stable vertebral fractures (lumbar spine) with vertebral body damage – relieving the anterior edges of vertebrae
- Moderate degeneration / muscular insufficiency of the spine
- Moderate lumbosciatica with muscular deficiencies in cases of intervertebral disk protrusions / herniated disks
- Facet syndrome
- Moderate (pseudo) radicular lumbar syndrome
- Lumbar spinal canal stenosis (non-surgical, post-operative)
- Spondylolysis / spondylolisthesis, grade II
- Condition after herniation (conservative / post-operative)
- Condition following discectomy
- Condition following spinal fusion / kyphoplasty (oneor more levels)
- Osteoporosis

BENEFITS AND MODES OF ACTION

- Integrated circular corset stays ensure stabilization and straightening
- Abdominal pad for increased compression of the abdomen to straighten the spine
- Supports natural lumbar lordosis using an anatomically shaped back pad

THE ADVANTAGES OF BAUERFEIND

- + Multi-stage treatment
- + Gentle on the skin and breathable

ORDERING OPTIONS

Colors: titan

Accessories / spare parts: see page 112 and price list

SIZE	1	2	3	4	5	6	ORDER NO.
Circumference trunk (cm)	70 - 80	80-90	90-105	105-115	115-125	125-140	122 8563 1000 00°
° = Size (1-6)							

Spinova[®] Unload Classic

Moderate lordosis reduction and relief of the lumbar spine using removable bridging frame









Integrated circular corset stays, can be adjusted

Removable bridging frame

Lordosis reduction

INDICATIONS

- Vertebral fractures (lumbar spine) with vertebral body damage, relieving the posterior edges of vertebrae
- Spondylolysis / spondylolisthesis, grades I–II
- Moderate lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Lumbar spinal canal stenosis, also with paresis (non-surgical, post-operative)
- Moderate (pseudo) radicular lumbar syndrome
- Facet syndrome
- Condition following discectomy

BENEFITS AND MODES OF ACTION

- Integrated circular corset stays ensure stabilization and straightening
- Corrects lumbar lordosis using a selfadjusting bridging frame
- Relieves pressure on the posterior edges of the lumbar vertebrae

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, long

- short: body height < 170 cm
- long: body height > 170 cm

Accessories/spare parts:

see page 112 and price list

SPINE 8 NECK

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-105	105–115	115–125	125–140	122 8562 100+ 00°	


+ = Design (1) short, (2) long ° = Size (1-6)

experts.bauerfeind.com/spinova-unload-classic

Spinova[®] Unload Plus

Lordosis reduction and increased relief of the lumbar spine with removable bridging frame







Lordosis reduction

Integrated circular corset stays, can be adjusted

Removable bridging frame

Abdominal pad

INDICATIONS

- Vertebral fractures (lumbar spine) with vertebral body damage relieving the posterior edges of vertebrae
- Spondylolysis/spondylolisthesis, grades I-II
- Most severe lumbar deformity in cases of facet syndrome / osteoarthritis
- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Lumbar spinal canal stenosis, also with paresis (non-surgical, post-operative)
- Most severe (pseudo) radicular lumbar syndrome
- Foraminal stenosis

BENEFITS AND MODES OF ACTION

- Integrated circular corset stays ensure stabilization and straightening
- Corrects lumbar lordosis using a selfadjusting bridging frame
- Relieves pressure on the posterior edges of the lumbar vertebrae
- Greater straightening thanks to increased abdominal pressure exerted by the abdominal pad

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, long

- short: body height < 170 cm
- long: body height > 170 cm

Accessories / spare parts:

see page 112 and price list

SIZE	1	2	3	4	5	6	7	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-100	100-110	110-120	120-130	130-140	122 8543 100+ 00°	-

+ = Design (1) short, (2) long ° = Size (1 - 7)

experts.bauerfeind.com/spinova-unload-plus

Spinova[®] Immo Plus

Immobilization of the lumbar spine thanks to mobilization function using a removable shell







Integrated circular corset stays, can be adjusted

Individually adjustable shell system in two sizes

INDICATIONS

- Vertebral fractures (lumbar spine) with significant vertebral body damage – relieving the anterior and / or posterior edges of vertebrae
- Tumors (metastases)
- Condition following spinal fusion (one or more levels)
- Most severe (pseudo) radicular lumbar syndrome
- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Lumbar spinal canal stenosis with paresis (non-surgical, postoperative)
- Foraminal stenosis
- Osteochondrosis
- Severe degeneration / muscular insufficiency
- Severe lumbar deformity in cases of facet syndrome/osteoarthritis
- Condition following herniation (non-surgical / post-operative)
- Condition following discectomy
- Spondylolysis/spondylolisthesis, grades I–III

Optional abdominal pad for increased compression of the abdomen

BENEFITS AND MODES OF ACTION

- The lumbar spine is stabilized, relieved and its movement limited using a shell surrounding the torso combined with a circular pelvic frame
- Individual adjustment and the perfect fit thanks to two-part shell system
- Together with the shell, the Spinova tensioning strap system forms a single functional unit; circumferential forces can be applied individually and to suit the indication
- Greater relief of the spine thanks to increased abdominal pressure exerted by the abdominal pad
- Integrated corset stays provide additional stabilization and straighten the spine



Immobilization

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, standard

- short: body height < 170 cm
- standard: body height > 170 cm

Accessories / spare parts:

see page 112 and price list

+ = Design (0) standard, (1) short ° = Size (1 - 7)

experts.bauerfeind.com/spinova-immo-plus

SPINE 8 NECK

Spinova[®] Immo Plus Classic

Immobilization of the lumbar spine thanks to mobilization function and removable shell





Individually adjustable shell system

Optional abdominal pad for increased compression of the abdomen

frame

spine

INDICATIONS

- Stable vertebral fractures (lumbar spine) with vertebral body damage
- Relieving the anterior and / or posterior edges of vertebrae
- Most severe lumbosciatica with muscular deficiencies in cases of intervertebral disk protrusions / herniated disks
- Lumbar spinal canal stenosis, also with paresis (non-surgical, post-operative)
- Tumors (metastases)
- Severe degeneration / muscular insufficiency of the spine
- Spondylolysis / spondylolisthesis, grades II–III
- Facet syndrome/spondylitis
- Most severe (pseudo) radicular lumbar syndrome
- Foraminal stenosis
- Condition following herniation (non-surgical/post-operative)
- Condition following discectomy
- Condition following spinal fusion / kyphoplasty (oneor more levels)
- Osteochondrosis



• The lumbar spine is stabilized, relieved and

its movement limited using a shell frame

· Individual adjustment and the perfect fit

thanks to two-part shell system

• Greater relief of the spine thanks to increased abdominal pressure exerted by

Integrated spring steel stays provide

additional stabilization and straighten the

the abdominal pad

for the torso combined with a circular pelvic

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Multi-stage treatment

ORDERING OPTIONS

Colors: titan

Designs: short, standard

- short: body height < 170 cm
- standard: body height > 170 cm

Accessories / spare parts: see page 112 and price list

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	70 - 80	80-90	90-105	105–115	115–125	125-140	122 8564 100+ 00°	-

+ = Design (0) standard, (1) short ° = Size (1-6)

experts.bauerfeind.com/spinova-immo-plus-classic

Spinova[®] Immo

Immobilization of the lumbar spine







Individually adjustable shell system

Immobilization

INDICATIONS

- Immobilization of the lumbar spine (pre-operative, post-operative, post-traumatic, inflammatory conditions)
- Non-surgical treatment of simple and stable fractures

BENEFITS AND MODES OF ACTION

- The lumbar spine is stabilized, relieved and its movement limited using a shell surrounding the torso combined with a circular pelvic frame
- Individual adjustment and the perfect fit thanks to two-part shell system
- Together with the shell, the Spinova tensioning strap system forms a single functional unit; circumferential forces can be applied individually and to suit the indication

5

6

4

THE ADVANTAGES OF BAUERFEIND

+ Digital measurement

ORDERING OPTIONS

Colors: titan

ORDER NO.

7

100-110 110-120 120-130 130-140 122 8540 1000 00°

Accessories / spare parts: see page 112 and price list

SPINE & NECK

1	

° = Size (1-7)

Circumference trunk (cm)

SIZE

1

70-80

2

80-90

3

90-100

SecuTec[®] Dorso

Lordosis reduction and relief of the lumbar spine using a bridging frame and corset stays







Removable bridging frame

Anatomically shaped corset stays



grid system

INDICATIONS

- Spondylolisthesis (with recurring lumbosciatica), grades I – III
- Condition following herniation (e.g. post-operative)
- Condition following discectomy
- Vertebral fractures (lumbar spine) with vertebral body damage
- Osteoporosis
- Moderate to severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Severe (pseudo) radicular lumbar syndrome
- Moderate to severe lumbar deformity in cases of facet syndrome / osteoarthritis
- Lumbar spinal canal stenosis

BENEFITS AND MODES OF ACTION

- Relieves and straightens the lumbar spine thanks to the combination of corset stays, anatomically contoured knitted fabric and tensioning straps
- Bridging frame ensures lordosis reduction

THE ADVANTAGES OF BAUERFEIND

+ Easy to handle

ORDERING OPTIONS

Colors: titan

SIZE	1	2	3	4	5	6	ORDER NO.	
Circumference trunk (cm)	75-85	85-95	95–105	105–115	115–125	125–135	121 8474 0080 00°	

° = Size (1-6)

SecuTec[®] Lumbo

Lordosis reduction and relief of the lumbar spine using a bridging frame







Removable bridging frame

Abdominal pad

INDICATIONS

- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Most severe (pseudo) radicular lumbar syndrome
- Spondylolysis/spondylolisthesis, grades I-III
- Moderate to severe lumbar deformity in cases of facet syndrome/osteoarthritis
- Lumbar spinal canal stenosis
- Moderate to severe degeneration / muscular insufficiency
- Vertebral fractures (lumbar spine) with vertebral body damage
- Condition following herniation, e.g. post-operative
- Condition following discectomy

BENEFITS AND MODES OF ACTION

- Lordosis-reducing effect thanks to the combination of removable bridging frame, anatomically contoured knitted fabric and tensioning straps
- Greater straightening thanks to increased abdominal pressure exerted by the abdominal pad

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Digital measurement

ORDERING OPTIONS

Colors: titan

ORDER NO.

121 8460 108+ 00°

Designs: short, long

- short: body height < 170 cm
- long: body height > 170 cm

Accessories / spare parts: see page 111 and price list

c	

Circumference trunk (cm) + = Design (1) short, (2) long ° = Size (1-6)

SIZE

experts.bauerfeind.com/secutec-lumbo

2

80-90

3

90-100

4

5

100-110 110-120 120-130

6

1

70-80

SofTec[®] Lumbo

Immobilization of the lumbar spine thanks to mobilization function using removable shell







Short version of half-shell (29 cm)

Long version of half-shell (35 cm)

INDICATIONS

- Spondylolysis / spondylolisthesis, grades II – III
- Most severe lumbosciatica with muscular deficiencies in cases of significant intervertebral disk protrusions / herniated disks
- Most severe (pseudo) radicular lumbar syndrome
- Degenerative facet syndrome
- Condition following herniation, e.g. post-operative
- Condition following discectomy
- Condition following spinal fusion / kyphoplasty (one or more levels)
- Lumbar spinal canal stenosis
- Vertebral fractures (lumbar spine) with vertebral body damage
- Tumors (metastases)
- Severe degeneration / muscular insufficiency



LPT strap system that requires minimal effort

BENEFITS AND MODES OF ACTION

- Stabilization of the lumbar spine: different degrees of stabilization can be achieved with the targeted combination of different orthosis modules (knitted fabric, two grades of CFRP stays, half-shell)
- Convenient for the patient to put on with minimal effort
- High level of wearing comfort
- Greater relief of the spine thanks to increased abdominal pressure exerted by the abdominal pad



Abdominal pad

THE ADVANTAGES OF BAUERFEIND

- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, long, two heights

- Dorsal height of the plastic half-shell: short 29 cm, long 35 cm
- Stabilization: short L5-S1, long L3-S1
- short: body height < 170 cm
- long: body height > 170 cm

Accessories / spare parts: see page 111 and price list

SIZE	1	2	3	4	HALF-SHELL	ORDER NO.	
Circumforonce hin (om)			105-120		short	122 8480 0081 00°	\sim
Circumference hip (cm)	75 – 90	90–105	105–120	120–135	long	122 8480 0082 00°	

• = Size (1-4)

experts.bauerfeind.com/softec-lumbo

DorsoTrain®

Active straightening and support of the spine







Reclinator made from fiberglass-reinforced plastic

Non-elastic functional zones

INDICATIONS

- Mild degeneration / muscular insufficiency
- Osteopenia, osteoporosis
- Stable, osteoporotic vertebral fractures

BENEFITS AND MODES OF ACTION

- Straightening and support of the spine using a semi-flexible reclinator combined with non-elastic functional zones of the bodysuit
- The anatomically shaped reclinator stabilizes the spine and prevents the development of a hunched back
- Individual adjustment of the degree of stabilization using tensioning straps
- Easy and quick to put on using zipper at the front and comfortable, high gusset fastening

THE ADVANTAGES OF BAUERFEIND

+ Gentle on the skin and breathable

ORDERING OPTIONS

Colors: ice-blue

Designs: short, long

- short: body height < 170 cm
- long: body height > 170 cm

Accessories / spare parts: see page 111 and price list

SPINE & NECK

SIZE	1	2	3	4	ORDER NO.
Making-up	S	М	L	XL	110 8455 008+ 00°
+ = Design (1) short, (2)	long				

° = Size (1 – 4)

Spinova[®] Osteo

Stabilization and active straightening of the spine, with free abdomen and chest area













Free chest and abdominal area

Pelvic support

Reclinator

Strap system

Active principle

INDICATIONS

Active straightening to relieve and correct the lumbar and thoracic spine as well as to restrict movement in the sagittal plane, e.g. in the case of

- Stable, osteoporotic vertebral fractures
- Osteoporosis, osteopenia
- Painful hunched back
- Muscular insufficiency

BENEFITS AND MODES OF ACTION

- Increased sense of stability thanks to padded crossed straps at the sides with a reminder effect
- Straightening effect can be regulated using a special strap system which allows individual pressure distribution
- Straightens using an aluminum reclinator combined with tensioning straps and a pelvic support
- Free abdomen and chest area as well as good positional stability of the orthosis due to low-position pelvic support
- Adaptable strap system thanks to heightadjustable dorsal strap fastenings
- Donning with minimal effort using the hand straps and finger loops

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Tool-free fitting

ORDERING OPTIONS

Colors: titan

Optional:

If the measurements are between two sizes, choose a smaller reclinator/bigger support.

Five standard sizes: S2, S3, M3, M4, L4

18 size combinations of reclinator (S-L) and support (1-6)

Accessories / spare parts:

see page 112 and price list

Spinova Osteo

SIZE	S 2	S 3	М3	M4	L4	ORDER NO.	7
Circumference back (cm)	50-58	50-58	59–67	59–67	68-76		
Circumference pelvis (cm)	80-90	90-100	90-100	100–110	100-110	122 8568 0000 0°°	90° 1
°° number	12	13	23	24	34		

Spinova Osteo with partially assembled reclinator

SIZE	S	М	L	ORDER NO.	
Length back (cm)	50-58	59-67	68 - 76	122 8570 0000 00°	
° = Size (1) S, (2) M, (3) L		Please no	ote: Choose a smaller reclinator if you are between sizes.		_

Spinova Osteo with partially assembled pelvic support

SIZE	1	2	3	4	5	6	ORDER NO.
Circumference pelvis (cm)	70-80	80-90	90-100	100-110	110-120	120-130	122 8569 0000 00°
° = Size (1-6)		Please no	ote: Choose a	larger pelvi	c support if y	ou are between si	zes.

experts.bauerfeind.com/spinova-osteo

SPINE 8 NECK


83

SofTec[®] Dorso

Increased stabilization of the thoracic and lumbar spine









Shoulder straps

Aluminum reclinator

Integrated cross

INDICATIONS

- Non-surgical treatment of simple and stable fractures (up to TH 8)
- Condition following discectomy
- Condition following spinal fusion
- Osteoporosis
- Tumors (metastases) or inflammatory processes (spondylitis)

BENEFITS AND MODES OF ACTION

- Straightens the spine and stabilizes the upper body using an adjustable aluminum reclinator
- Protects against damaging, incorrect movements
- Pressure distribution and straightening can be regulated using tensioning straps
- A stabilizing cross made of aluminum that can be fixed to the reclinator increases the straightening effect of the orthosis and provides additional restriction in the end range of motion (extension / flexion; bending to the side)

THE ADVANTAGES OF BAUERFEIND

+ Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: short, long

- short: body height < 170 cm
- long: body height > 170 cm

Accessories / spare parts: see page 111 and price list

PROVEN IN A STUDY^Δ

• Provides relief, pain reduction and mobilization

SIZE	1	2	3	4	ORDER NO.	
Circumference hip (cm)	75-90	90-105	105-120	120–135	122 8490 008+ 00°	
+ = Design (1) short, (2) lor	ng					

° = Size (1 – 4)

^a Bauerfeind, non-interventional study, unpublished data

experts.bauerfeind.com/softec-dorso

CerviLoc[®] | CerviLoc[®] S

Stabilization and relief of the cervical spine





Side view

INDICATIONS

- Pain in the cervical spine (cervical syndrome)
- Whiplash

CerviLoc S

- Severe pain in the cervical spine (severe cervical syndrome)
- Moderate whiplash
- Blockage

BENEFITS AND MODES OF ACTION

- Targeted relief of the cervical spine
- The integrated flexible plastic reinforcement of the CerviLoc S provides increased support
- Comfortable to wear thanks to soft, anatomically shaped foam with a cut-out for the chin

APPLICATION INSTRUCTIONS

• Three different heights, each with two different circumferences, take into account individually varying neck shapes

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: beige

85

SIZE	1	2	3	4	5	6	ORDER NO.	
Height (cm)	6.5	8.5	10	6.5	8.5	10	120.050, 0010.00%	
Circumference neck (cm)	32-38	34-40	34-40	36-44	36-44	38-48	120 959+ 0010 00°	آدب

+ = Design (0) CerviLoc, (1) CerviLoc S ° = Size (1-6)

experts.bauerfeind.com/cerviloc | experts.bauerfeind.com/cerviloc-s

FOR TARGETED HAND, ELBOW AND SHOULDER TREATMENT



INDICATION GUIDE FOR ARM AND SHOULDER TREATMENT

PRODUCT		ManuTrain°	EpiTrain*	OmoTrain*	OmoTrain [®] S	ManuLoc*	ManuLoc° long
							- AN
PRODUCT FEATURI	E	Relief and stabilization of the wrist	Relief and stabilization of the elbow joint	Relief and stabilization of the shoulder joint	Relief and stabilization of the shoulder joint	Immobilization of the wrist	Immobilization of the wrist and parts of the forearm
PAGE		90	98	100	101	91	92
NDICATIONS							
Irritation (post-traumatic / post-ope	erative)	•	•	•	•	•	•
Pre-operative care						•	٠
(Relapse) prevention		٠	•	•	•		
Sprains		•				•	•
Contusions						•	•
nstability, e.g. radioulnar	subluxation	•					
	mild	•	•	•	•		
Osteoarthritis / arthritis	moderate					•	•
	severe						
Collateral ligament injurie at the thumb, e.g. skier's	es thumb						
Ligament injuries in the fi hand area	inger /						
Paralysis							
	mild	٠					
Tendovaginitis tenosynovitis)	moderate					٠	٠
	severe						
	mild					•	
Carpal tunnel syndrome	moderate						•
	severe						
Tendomyopathy			•				
Lateral /medial humeral epicondylitis			•				
Rotator cuff rupture (complex)							
Subcapital humerus fract	ure						
SLAP lesion							
Shoulder replacement su	rgery						
Shoulder dislocation							
Acromioclavicular joint in	iury (AC rupture)						

		STABILIZATION	I				RELIEF
ManuLoc° long Plus	RhizoLoc*	ManuLoc° Rhizo	ManuLoc® Rhizo long	ManuLoc® Rhizo long Plus	OmoLoc°	SecuTec [®] Omo	EpiPoint [®]
- AN	D.	13	and a	- And	and a	AU	-
mmobilization of the vrist, parts of the orearm and the ingers	Stabilization of the metacarpophalangeal and carpometacarpal joints of the thumb	Immobilization of the wrist, with thumb brace	Immobilization of the wrist and parts of the forearm, with thumb brace	Immobilization of the entire hand and parts of the forearm	Immobilization of the shoulder	Immobilization of the shoulder with individual arm positioning and mobilization function	Targeted relief of the tendon insertions at the elbow
93	94	95	96	97	102	103	99
•	•	•	•	•	•	•	
•		acute	•	•	•	•	
•	•	•	•	•			٠
•	•	•	•	•			
	•						
		•	•				
	•	•	•	•			
•		•		•			
•		•	•	•			
		•					
•			•				
							•
						•	
						•	
						•	
					•	•	

ManuTrain®

Relief and stabilization of the wrist





Two anatomically contoured pads

INDICATIONS

- Chronic, post-traumatic or post-operative irritation
- Distortions (wrist injuries, e.g. sprains)
- Osteoarthritis (joint wear)
- Instability, e.g. radioulnar subluxation
- Tendovaginitis (tenosynovitis)
- (Relapse) prevention

BENEFITS AND MODES OF ACTION

- Stabilization and activation of the wrist through even, circular compression
- Viscoelastic, anatomically contoured pads relieve pressure on nerves and vessels
- Indication-specific stabilization of the hand thanks to shapeable inner hand stay
- Inner hand stay and Velcro fastening provide greater relief to the wrist
- The seamless knitted fabric between the thumb and index finger enables gripping of objects without irritation

APPLICATION INSTRUCTIONS

• The shapeable inner hand stay allows functional adjustment of the hand for various indications

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement

ORDERING OPTIONS

Colors: beige, black, titan **Designs:** right, left

Accessories / spare parts: see page 114 and price list

SCOPE OF DELIVERY

Including donning aid

PROVEN IN A STUDY^Δ

• Significant alleviation of problems after 2 weeks

SIZE	1	2	3	4	5	6	ORDER NO.
Circumference wrist (cm)	14–15	15–16	16–17	17–18	18–19	19-20	110 5150 30*0 +0°

* = Color (1) beige, (7) black, (8) titan

+ = Design (6) right, (7) left ° = Size (1-6)

^a Spallek M., Akute Handgelenk- und Unterarmbeschwerden; Schnell wieder fit mit muskelaktivierender Bandage; Ärztliche Praxis, 44. Jahrgang, Nr. 20, S. 12–14. 1992

experts.bauerfeind.com/manutrain

ManuLoc[®]

Immobilization of the wrist







Metacarpal stay and stays on the sides

Edges finished using ultrasonic welding



Velcro fastening with retraction stopping device

INDICATIONS

- Post-traumatic or post-operative irritation
- Pre-operative care
- Carpal tunnel syndrome (mild)
- Tenosynovitis (moderate)
- Sprains/contusions
- Osteoarthritis (moderate)

BENEFITS AND MODES OF ACTION

- Securing the wrist during both extension/flexion and radial/ulnar movements thanks to aluminum stays
- The orthosis prevents incorrect positioning of the distal and proximal parts of the wrist
- Anatomically shaped metacarpal stay stabilizes the wrist in a functional position and can be individually adjusted
- The design of the orthosis enables functional gripping movements, thus maintaining mobility
- Practical Velcro fastenings and retraction stopping devices on the straps make it easy for the patient to put on the orthosis using only one hand
- Great freedom of movement for the thumb thanks to large thumb openings on both sides - the orthosis can therefore be worn on the left and the right
- Edges finished using ultrasonic welding for increased durability and an improved look

APPLICATION INSTRUCTIONS

• Metacarpal stay with ventilation slots and stays at the sides made of aluminum are individually shapeable

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: universal (thumb openings on both sides)

SIZE	0	1	2	3	ORDER NO.	. ///.
Circumference wrist (cm)	12–15	14–18	17-21	20-26	120 5340 2080 00°	



• = Size (0-3)

ManuLoc[®] long

Immobilization of the wrist and parts of the forearm





• Post-traumatic or post-operative

irritation, e.g. following treatment

with a plaster cast in the case of

stays on the sides

INDICATIONS

radius fractures

• Pre-operative care

(moderate)

Tenosynovitis

• Carpal tunnel syndrome

resting and at night)

Osteoarthritis (moderate)

• Sprains / contusions

(moderate, with pain while





welding



Velcro fastening with retraction stopping device

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: universal (thumb openings on both sides)

Accessories / spare parts: see page 114 and price list

22% improvement in functional	-
ity during everyday activities s	ix

PROVEN IN A STUDY^Δ

ies six weeks after surgery



Edges finished using ultrasonic

BENEFITS AND MODES OF ACTION

- Securing the wrist during both extension / flexion and radial / ulnar movements thanks to extended aluminum stays
- Limitation of rotational movements in the wrist thanks to the orthosis' longer lever
- the distal and proximal parts of the wrist
- thanks to large thumb openings on both sides - the orthosis can therefore be worn on the left and the right
- whilst maintaining functional gripping movements thanks to shapeable metacarpal stay
- Practical Velcro fastenings and retraction for the patient to put on the orthosis using only one hand
- Edges finished using ultrasonic welding for increased durability and an improved look; sensitive skin in the elbow area is softly cushioned / enclosed

APPLICATION INSTRUCTIONS

• The extended stays on the sides and the metacarpal stay with ventilation slots made of aluminum are individually shapeable

(a) 21 cm

SIZE ORDER NO. Π 2 3 1 120 5340 2082 00° Circumference wrist (cm) 12 - 1514 - 1817-21 20-26

• = Size (0-3)

^a Dr. med. Peter Katzmaier, POSTOPERATIVE THERAPIE DER DISTALEN RADIUSFRAKTUR MIT DER ORTHESE ManuLoc long Bauerfeind Life 3-2016. The full article can be accessed at bauerfeind-life.de

experts.bauerfeind.com/manuloc-long

- The orthosis prevents incorrect positioning of
- Great freedom of movement for the thumb
- Immobilization in a functional position
 - stopping devices on the straps make it easy



ManuLoc[®] long Plus

Immobilization of the wrist. parts of the forearm and the fingers







Shapeable hand support

Edges finished using ultrasonic welding

INDICATIONS

stays on the sides

- Post-traumatic care (multiple, complex injuries in the hand / finger area)
- Pre-operative and post-operative, e.g. Dupuytren's contracture, metacarpal fractures
- Carpal tunnel syndrome (severe - support at night)
- Tenosynovitis (severe)
- Sprains/contusions
- Ligament injuries in the finger/hand area

- Stabilizes the wrist, forearm and fingers, while the thumb remains freely moveable
- Adjustable rest positions in Intrinsic Plus position or normal position to reduce pain and prevent damage associated with immobilization
- Securing the wrist during both extension/flexion and radial/ulnar movements thanks to extended aluminum stays
- The wide finger strap relieves the ligaments and tendons while the fingers are held in position.
- Great freedom of movement for the thumb thanks to large thumb openings on both sides - the orthosis can therefore be worn on the left and the right
- Practical Velcro fastenings and retraction stopping devices on the straps make it easy for the patient to put on the orthosis using only one hand
- Immobilization in a functional position thanks to shapeable metacarpal stay
- Edges finished using ultrasonic welding for increased durability and an improved look; sensitive skin in the elbow area is softly cushioned / enclosed

APPLICATION INSTRUCTIONS

• The extended stays on the sides and the metacarpal stay with ventilation slots made of aluminum are individually shapeable

SIZE	0	1	2	3	ORDER NO.
Circumference wrist (cm)	12-15	14–18	17-21	20-26	128 5341 0080 00°

• = Size (0 – 3)



Velcro fastening with retraction stopping device

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: universal (thumb openings on both sides) HAND, ELBOW & SHOULDER

Accessories / spare parts:

see page 114 and price list

93





RhizoLoc®

Stabilization of the metacarpophalangeal and carpometacarpal joints of the thumb









Velcro fastening with retraction stopping device

Velcro tab for mobilization Aluminum frame

INDICATIONS

- Post-operative and post-traumatic irritation
- Collateral ligament injuries, e.g. skier's thumb
- Osteoarthritis of the carpometacarpal joint of the thumb (mild)
- Sprains/contusions

BENEFITS AND MODES OF ACTION

- Stabilization of the metacarpophalangeal and carpometacarpal joints of the thumb
- Excellent stabilization of the metacarpophalangeal joint of the thumb thanks to aluminum frame
- Restrictable release of movement at the metacarpophalangeal joint of the thumb (dynamization) thanks to Velcro tab
- Practical Velcro fastenings and retraction stopping devices on the straps make it easy for the patient to put on the orthosis using only one hand
- High degree of wearing comfort based on breathable materials with air holes
- Freedom of movement retained in the wrist thanks to the slim strap and small deflection ring

APPLICATION INSTRUCTIONS

• Individually adaptable aluminum frame

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting

ORDERING OPTIONS

Colors: titan **Designs:** right, left

PROVEN IN A STUDY[∆]

• Pain reduction by immobilizing the carpometacarpal joint of the thumb

(up to MCP joint I) in cases of polyosteoarthritis in the fingers

SIZE	1	2	ORDER NO.	
Circumference wrist (cm)	12.5–16	16-20	120 5220 1000 - 00	
Circumference thumb (cm)	5 – 7	6.5-9	120 5330 1080 +0°	21

+ = Design (6) right, (7) left ° = Size (1-2)

^a Wanivenhaus A., Die Behandlung der Fingerpolyarthrosen; Rheumatologie und Orthopädie, 1, 10–12, 2013

experts.bauerfeind.com/rhizoloc

ManuLoc[®] Rhizo

Immobilization of the wrist. with thumb brace











Edges finished using ultrasonic welding

Velcro fastening with retraction stopping device

INDICATIONS

stays on the sides

- · Acute, post-traumatic, and post-operative irritation in the wrist and thumb ray
- Pre-operative care
- Tenosynovitis at the radial side (moderate)
- Osteoarthritis of the carpometacarpal joint of the thumb (severe)
- Carpal tunnel syndrome (mild)
- Sprains / contusions
- Ligament injuries in the finger / hand area



BENEFITS AND MODES OF ACTION

- The orthosis prevents incorrect movements of the proximal and distal parts of the wrist as well as of the metacarpophalangeal and carpometacarpal joints of the thumb.
- Immobilization of the wrist thanks to shapeable aluminum stays
- Stabilization of the wrist in a functional position thanks to anatomically shaped metacarpal stay
- Practical Velcro fastenings and retraction stopping devices on the straps make it easy for the patient to put on the orthosis using only one hand.
- Edges finished using ultrasonic welding for increased durability and an improved look

APPLICATION INSTRUCTIONS

- The metacarpal stay with ventilation slots, the stays at the sides and the frame made of aluminum are individually shapeable.
- Over the course of the treatment, the thumb can be increasingly allowed movement and the thumb brace can eventually be removed completely.

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 114 and price list

SIZE ORDER NO. 0 2 3 20 - 26120 5345 0080 +0° Circumference wrist (cm) 12-15 14-18 17-21

+ = Design (6) right, (7) left ° = Size (0-3)



ManuLoc[®] Rhizo long

Immobilization of the wrist and parts of the forearm, with thumb brace





Metacarpal stay and stays at the sides

Individually adjustable thumb immobilize

INDICATIONS

- Post-traumatic or post-operative irritation, e.g. following treatment with a plaster cast in the case of scaphoid fractures
- Pre-operative care
- Tenosynovitis at the radial side (moderate)
- Osteoarthritis of the carpometacarpal joint of the thumb (severe)
- Sprains / contusions
- Carpal tunnel syndrome (moderate)
- Collateral ligament injuries, e.g. skier's thumb



Aluminum frame

BENEFITS AND MODES OF ACTION

- The orthosis prevents incorrect movements of the proximal and distal parts of the wrist as well as of the metacarpophalangeal and carpometacarpal joints of the thumb
- Limitation of rotational movements in the wrist thanks to the orthosis' longer lever
- Immobilization of the wrist thanks to shapeable aluminum stays
- Stabilization of the wrist in a functional position thanks to anatomically shaped metacarpal stay
- Practical Velcro fastenings and retraction stops on the straps make it easy for the patient to put on the orthosis using only one hand
- Edges finished using ultrasonic welding for increased durability and an improved look; sensitive skin in the elbow area is softly cushioned / enclosed

APPLICATION INSTRUCTIONS

- The metacarpal stay with ventilation slots, the stays at the sides and the frame made of aluminum are individually shapeable
- Over the course of the treatment, the thumb can be increasingly allowed movement and the thumb brace can eventually be removed completely



Edges finished using ultrasonic welding



Velcro fastening with retraction stopping device

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 114 and price list

SIZE	0	1	2	3	ORDER NO.		
Circumference wrist (cm)	12 – 15	14–18	17 – 21	20-26	120 5345 0082 +0°	2	

+ = Design (6) right, (7) left ° = Size (0-3)



 $\Delta > 21 \text{ cm}$

ManuLoc[®] Rhizo long Plus

Immobilization of the entire hand and parts of the forearm







Shapeable finger support

INDICATIONS

stays at the sides

- Post-traumatic care (multiple, complex injuries in the hand / finger area)
- Pre-operative and post-operative, e.g. Dupuytren's contracture, metacarpal fractures
- Complex ligament injuries in the finger / hand area due to trauma
- Paralysis
- Tenosynovitis (severe)
- Sprains / contusions



Individually adjustable thumb immohilizei

BENEFITS AND MODES OF ACTION

damage due to immobilization

including the thumb

aluminum stays

cushioned / enclosed

healing process

APPLICATION INSTRUCTIONS

hand

· Stabilizes the wrist, forearm and fingers

• Adjustable positions to reduce pain and prevent

• Immobilization of the wrist thanks to shapeable

stopping devices on the straps make it easy for

the patient to put on the orthosis using only one

• Practical Velcro fastenings and retraction

• Edges finished using ultrasonic welding for

sensitive skin in the elbow area is softly

• The finger support is supplied in a neutral position and can be fixed in the Intrinsic Plus position guickly and without the need for tools • The metacarpal stay with ventilation slots, the stays at the sides and the frame made of aluminum are individually shapeable

• The removable finger support and thumb brace aid physiotherapy and mobilization during the

increased durability and an improved look;



Edges finished using ultrasonic welding

Velcro fastening with retraction stopping device

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Comfort zone
- + Tool-free fitting
- + Multi-stage treatment
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Accessories / spare parts: see page 114 and price list

SIZE **ORDER NO.** 0 2 3 1 12 - 1517 - 21128 5347 0080 +0° Circumference wrist (cm) 14 - 1820 - 26

+ = Design (6) right, (7) left • = Size (0-3)

experts.bauerfeind.com/manuloc-rhizo-long-plus

A > 21 cm



97

EpiTrain[®]

Relief and stabilization of the elbow joint









Comfortable elbow zone



Epicon+ Pads

Epicon+ Pads on the product

Long upper arm knit

Train knit fabrio

INDICATIONS

- Post-traumatic or post-operative irritation, e.g. epicondylitis
- Osteoarthritis/arthritis (joint wear / joint inflammation)
- Tendinopathy (degenerative tendon condition)
- Tendomyopathy . (muscle-tendon condition)
- Prevention / relapse prevention

BENEFITS AND MODES OF ACTION

- The therapeutic compression of the specially textured surface of the high-low knitted fabric combined with the Epicon+ Pads provides noticeable relief and targeted pain reduction
- The viscoelastic Epicon+ Pads with epicondyle cut-outs redistribute the support's compression away from the epicondyles to the surrounding soft tissue and allow for a secure positioning of the support
- The soft and particularly breathable comfort zone in the elbow of the anatomically designed train knit allows for a high degree of wearing comfort
- The Epicon+ Pads relieve the epicondyles and the massage nubs with a firmer material component exert a soothing intermittent compression massage on the muscle and tendon attachment sites during movement

ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Comfort zone
- + Digital measurement

ORDERING OPTIONS

Colors: beige, black, titan

Versions: Identical versions for right and left

Optional: custom-made solutions in beige and titan, top band, velcro fastening

A 12

Accessories/spare parts:

see page 115 and price list

PROVEN IN A STUDY^A

• 50 % less pain after 14 days compared with the start of treatment using EpiTrain (control group: -19%)

SIZE	0	1	2	3	4	5	6	ORDER NO.
Largest forearm circumference (cm)	21-23	23-25	25-27	27–29	29-31	31-33	33-35	110 6160 40*0 00°
* = Color (1) beige, (7) black, (8) titan ° = Size (0-6)		made solutio m-made pro						

= Size (0 – 6)

^a Valle-Jones et al., Controlled trial of an ellbow support (EpiTrain) in patientes with acute painful conditions of the ellbow a pilot study, Curr. Med. Res. Opin., 12, S. 224–233, 1990

experts.bauerfeind.com/epitrain

EpiPoint[®]

Targeted relief of the tendon insertions at the elbow







Red warning section indicates excessive pressure

Five-point pad

INDICATIONS

- Lateral epicondylitis of the humerus (inflammation of the muscle origin, tennis elbow)
- Medial epicondylitis of the humerus (inflammation of the muscle origin, golfer's and thrower's elbow)
- (Relapse) prevention

BENEFITS AND MODES OF ACTION

- Targeted relief for areas subjected to excessive strain thanks to regulated compression on the forearm
- The viscoelastic five-point pad stimulates trigger points and enables highly targeted compression
- An elastic adjustable band allows individual continuous pressure adjustment
- The adjustable band's red warning section indicates excessive pressure

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: black, titan **Designs:** universal

SCOPE OF DELIVERY

Also in a container (12x)

PROVEN IN A STUDY^Δ

• 23% less pain (in patients with positive extensor grip test)

SIZE	UNIVERSAL	ORDER NO.	
Circumference forearm (cm)	22.5-33	120 6350 10*0 000	

* = Color (7) black, (8) titan

^a Doctoral thesis: Struijs et al.; Frederik Hendrikplantsoen 74-2; NL-1052 XW Amsterdam

experts.bauerfeind.com/epipoint

OmoTrain[®]

Relief and stabilization of the shoulder joint





Removable Delta pad

INDICATIONS

- Chronic, post-traumatic or post-operative irritation
- (Relapse)prevention
- Osteoarthritis of the shoulder

BENEFITS AND MODES OF ACTION

- The compressive knitted fabric and the tensile force exerted by the strap system help center the glenohumeral joint
- The viscoelastic Delta pad with frictional nubs massages the shoulder's soft tissue during movement; joint perception is improved, the absorption of edemas and effusions is accelerated and pain is reduced
- The strap system leaves the chest area free
- The Delta pad can be optionally inserted with Velcro tabs and positioned as required

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: Identical design for right and left

Accessories / spare parts:

see page 115 and price list

PROVEN IN A STUDY^Δ

• Supports functional treatment in cases of rotator cuff lesions or tendinosis

SIZE	0	1	2	3	4	5	6	ORDER NO.
Circumference (cm)	20-22	22-24	24-26	26-29	29-32	32-36	36-40	110 7170 2080 00°
$^{\circ}$ = Size (0 - 6)								



OmoTrain[®] S

Relief and stabilization of the shoulder joint





Removable Delta pad

INDICATIONS

- Chronic, post-traumatic or post-operative irritation
- (Relapse) prevention
- Osteoarthritis of the shoulder

BENEFITS AND MODES OF ACTION

- The compressive knitted fabric and the tensile force exerted by the strap system improve the coadaptation of the humeral head to the glenoid cavity
- The viscoelastic Delta pad with frictional nubs massages the shoulder's soft tissue during movement; this relieves pain, improves joint perception and activates the joint-stabilizing muscles
- The Delta pad can be optionally inserted with Velcro tabs and positioned as required
- The direct strap system makes it quick and easy to put on the support

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Train knit fabric
- + Digital measurement

ORDERING OPTIONS

Colors: titan **Designs:** right, left

SIZE	0	1	2	3	4	5	6	ORDER NO.
Circumference (cm)	20-22	22-24	24-26	26-29	29-32	32-36	36-40	110 7171 0080 +0°

+ = Design (6) right, (7) left ° = Size (0-6)

OmoLoc[®]

Immobilization of the shoulder







Fastening with finger pocket



Wide hand / arm loops provide extensive support for the arm.

INDICATIONS

Soft, padded tabs

- Immobilization, e.g. after shoulder surgery
- Post-traumatic, pre-operative care
- Shoulder dislocation



With shoulder strap

- Immobilization of the shoulder joint in a 0°-degree abduction position (rotated inward), e.g. to maintain the surgical outcome
- The forearm is fixed at the required flexion angle of 90 100° and can therefore relieve pain
- Three arm loops that can be positioned individually provide a secure and perfect fit
- Secure adjustment according to different situations using the abdominal fastening with tensioning strap system
- More hygienic thanks to the design of the orthosis
- No shoulder strap required; this avoids incorrect posture and the related tension that occurs in the healthy shoulder; the shoulder strap is still available as a spare part

APPLICATION INSTRUCTIONS

• Fixing of the forearm in the correct position for medical treatment – 90–100°

THE ADVANTAGES OF BAUERFEIND

- + Easy to handle
- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Optional: Shoulder strap available as spare part

Accessories / spare parts: see page 115 and price list

SIZE	1	2	ORDER NO.
Circumference trunk (cm)	up to 100	100–130	120 7360 1120 +0°
+ = Design (6) right, (7) left			

° = Size (1 – 2)

experts.bauerfeind.com/omoloc

SecuTec[®] Omo

Immobilization of the shoulder with individual arm positioning and mobilization function







Abduction / adduction

Rotation

INDICATIONS

- Rotator cuff rupture (complex)
- Subcapital humerus fracture
- SLAP lesion
- Shoulder replacement surgery
- Shoulder dislocation
- Post-traumatic, post-operative irritation
- Pre-operative care
- Acromioclavicular joint injury (AC rupture)

BENEFITS AND MODES OF ACTION

- Immobilization of the shoulder joint with freely adjustable abduction position to relieve the tendons
- Pain reduction through restriction
- Mobilization by allowing inward or outward rotation at the push of a button
- The height of the arm support can be adjusted to the length of the patient's arm; the angle of inclination is adjusted automatically
- Relief of the hand using an extendible hand support
- Secure and individual adjustment using an abdominal fastening with tensioning strap system
- Enabling rotation provides relief and increases wearing comfort during the night. This should only be done upon consultation with a physician
- No shoulder strap required; this avoids incorrect posture and the related tension that occurs in the healthy shoulder; the shoulder strap is still available as a spare part

APPLICATION INSTRUCTIONS

- Can be adjusted in two planes
- The abduction angle can be individually adjusted up to 60°
- inward and outward rotation can be regulated at the touch of a button
- The angle of the pelvic frame can be adjusted to ensure that the orthosis fits the patient's body shape
- No shoulder strap required because the patient's arm can be supported on the pelvic frame

THE ADVANTAGES OF BAUERFEIND

- + Gentle on the skin and breathable
- + Digital measurement

ORDERING OPTIONS

Colors: titan

Designs: right, left

Optional: Shoulder strap available as spare part

Accessories / spare parts: see page 115 and price list

SIZE	UNIVERSAL	ORDER NO.	
Circumference trunk (cm)	70-140	121 7505 0120 +01	
+ = Design (6) right, (7) left			



ACCESSORIES AND SPARE PARTS

Spare parts are available within 3 to 4 working days.

AchilloTrain[®]

106

DES	CRIPTION	ITEM NO.		
01	Compensating cushion (for support sizes $1-2$)	129 1102 1000 001		
02	Compensating cushion (for support sizes 3–6)	129 1102 1000 003		
03	Heel cushion (for support sizes 1 – 2)	119 1103 0000 +01		
04	Heel cushion (for support sizes 3–6)	119 1103 0000 +02	01 02	03 04
			_	

+ = Design (6) right, (7) left

AirLoc®

DES	CRIPTION	ITEM NO.	
01	SofTec OA / AirLoc inflation bulb	129 4684 0070 000	<u>t</u>



MalleoLoc[®] L

DES	CRIPTION	SIZE	ITEM NO.
01	Strap set	universal	129 1328 1000 +00
02	Cushion	universal	129 1329 1000 +00
03	Pad		129 1332 1000 000

+ = Design (6) right, (7) left

MalleoLoc[®] L3

DES	CRIPTION	SIZE	ITEM NO.
01	Strap set	universal	129 1328 1000 +00
+ = D	lesign (6) right, (7) left		



MalleoLoc®

DES	CRIPTION	SIZE°	ITEM NO.
01	Strap, proximal	universal	129 1318 3123 +00
02	Strap, distal	universal	129 1318 3124 000
03	Buckle		129 1309 3120 000
04	Strap set without buckle cushion	universal	129 1302 3080 +00
05	Strap set with buckle cushion	1 – 2	129 1302 3080 +0°
06	Orthometer / goniometer		121 4674 0030 000



+ = Design (6) right, (7) left

ValguLoc®

DES	CRIPTION	SIZE°	ITEM NO.	
01	Strap set	1-3	129 1322 0000 +0°	
+ = [Design (6) right, (7) left			
				01

ValguLoc[®] II

DES	CRIPTION	SIZE	ITEM NO.
01	Strap set	universal	129 1322 0000 000

Donning aid

DESCRIPTION		ITEM NO.
01	Donning aid for supports	140 0222 0050 000
		01

TREATMENT AREA: KNEE, THIGH AND HIP

GenuTrain[®]

COLOR	ITEM NO.	
titan	119 4123 0000 000	



GenuTrain[®] S

DESCRIPTION		SIZE	ITEM NO.	
01	Strap set (for support sizes 1–2)	1	129 4132 30*0 001	
02	Strap set (for support sizes 3–4)	2	129 4132 30*0 002	
03	Strap set (for support sizes 5–7)	3	129 4132 30*0 003	
04	Joint stays	universal	129 4131 30*0 +00	



* = Color (1) beige, (7) black, (8) titan + = Design (6) right, (7) left



DESCRIPTION		SIZE	ITEM NO.			
01	Joint splints, lateral	universal	129 4602 1060 +00		4.5	000
02	Joint splints, medial	universal	129 4603 1060 +00	C		*****
03	Joint protection (sports protection)		129 4601 1070 000			00000
04	Limitation set (5 covering caps, 10 lock pins)		129 4607 0000 000	01 02	03	04

+ = Design (6) right, (7) left

SecuTec[®] Genu

DESCRIPTION		SIZE°	ITEM NO.
01	Strap set	1 – 7	129 4611 2070 00°
02	Frame cushion set	1 – 7	129 4613 1070 +0°
03	Strap cushion set		129 4617 1070 00"
04	Strap cushion, siliconized (x2)		129 4698 0000 00^
05	Limitation wedges (2 x frames + 2 x E0)		129 4616 1060 000
06	Frame	1 – 7	129 4650 1**# +0°
07	Joint pouch		129 4651 1070 000
08	Tibial edge cushion		129 4652 0070 000
09	Joint cover with quarter turn fastening (x1)		129 4645 1000 000
10	Quick-fit fastener set (x4)		129 4654 1000 000



GenuTrain[®] OA

DESCRIPTION		ITEM NO.		
01	Limitation set for sizes 1 and 2	129 4166 0070 001		
02	Limitation set for sizes 3 and 4	129 4166 0070 003		
03	Limitation set for size 5	129 4166 0070 005		



SofTec[®] OA

DES	SCRIPTION	SIZE	ITEM NO.		
01	Hinge – complete (with air cushion, pressure pad)	universal	129 4682 0070 +00		
02	Air cushion	universal	129 4683 0000 +00	-	
03	SofTec OA / AirLoc inflation bulb		129 4684 0070 000	9	
04	Limitation set (4 flexion stops, 3 extension stops)		129 8681 0060 000	01	02



+ = Design (6) right, (7) left

DES	CRIPTION	SIZE°	ITEM NO.
01	Strap set	1 – 7	129 4693 0000 00°
02	Frame cushion set (incl. gel cushion)	1 – 7	129 4694 0000 +0°
03	Strap cushion set	1 – 7	129 4695 0000 00°
04	Strap cushion, siliconized (x2)		129 4698 0000 00^
05	Frame (blue)	1 – 7	129 4696 005# +0°
06	OA joint pouch		129 4697 0000 000
07	Quick-fit fastener set (x4)		129 4654 1000 000
08	Limitation wedges (2 x frames + 2 x E0)		129 4616 1060 000
09	Allen key 3 mm		190 0007 0010 102



+ = Design (6) right, (7) left # = Design (0) distal, (3) proximal ^ = (1) for sizes 1 - 2, (3) for sizes 3 - 4, (5) for sizes 5 - 6, (7) for size 7

MOS-Genu

DES	CRIPTION	SIZE°	ITEM NO.
ADI	DITIONAL PARTS		
01	Popliteal cross-strap (for sizes 1–2)		121 4662 4050 001
02	Popliteal cross-strap (for size 3)		121 4662 4050 003
03	Popliteal cross-strap (for sizes 4–5)		121 4662 4050 004
04	Popliteal pads	1 – 5	121 4661 4090 00°
05	Joint ring	1 – 5	121 4664 4050 +0°
06	Spare parts set (screws)		121 4671 3000 000
07	Covering disks (x10)		121 4673 0000 000
08	Orthometer / goniometer		121 4674 0030 000
SHO	DRT VERSION		
09	Frame	1 – 5	121 4664 405# +0°
10	Cushion set	1 – 5	121 4665 4090 +0°
11	Strap set	1 – 5	121 4667 4050 00°
LON	IG VERSION		
12	Extension frame	universal	121 4669 405# +00
13	Strap for proximal extension frame, for sizes 1–3	1-3	121 4668 4053 00°
14	Strap for proximal extension frame, for sizes 4–5		121 4668 4053 004
15	Strap for distal extension frame, for sizes 1–3		121 4668 4054 001
16	Strap for distal extension frame, for sizes 4–5		121 4668 4054 004
17	Cushion for extension frame	universal	121 4666 409# +00

P	11	Ş	0000 0000 * h 4 k
01 02 03	04	05	06
		A	
07	08	09	10
	A	/	
11	12	13 14 15 16	17


+ = Design (6) right, (7) left # = Design (3) proximal, (4) distal

CoxaTrain[®]

DESCRIPTION

01	Covering cap	129 2583 1000 000
02	Gluteal pad	129 2585 1000 000
03	Trochanter pads / cushion	129 2586 1000 000
04	Joint limiting pins (2×)	129 2587 1000 000
05	Hinge with caps, short	129 2588 1001 000
06	Hinge with caps, long	129 2588 1002 000
07	Pads (2-part)	119 8416 0000 000

ITEM NO.



$\mathbf{SofTec}^{\circ} \mathbf{Coxa}$

DES	CRIPTION	SIZE°	ITEM NO.
01	Bag of small components		129 2521 0000 000
02	Thigh support, replacement part	universal	129 2521 0080 000
03	Pelvic knitted fabric	1 – 3	129 2521 0080 00°
04	Hinge, straight	universal	122 2523 1080 000
05	Hinge, offset	universal	122 2522 1080 001
06	Limitation set (5 covering caps, 10 lock pins)		129 2534 0000 000
07	Hip cup with assembly set	1 – 3	129 2535 0000 +0°
08	Allen key 3 mm		190 0007 0010 102



+ = Design (6) right, (7) left

PLEASE NOTE: A complete SofTec Coxa cannot be made from these spare parts

TREATMENT AREA: SPINE AND NECK

Fastening extension

DESCRIPTION	SIZE	ITEM NO.	
Back product fastening extension, short Can be used for the following product: SofTec Lumbo short	universal	129 8488 0082 000	
Back product fastening extension, long Can be used for the following products: SofTec Lumbo long, SofTec Dorso long + short, SecuTec Dorso, SecuTec Lumbo, LumboLoc Forte, LordoLoc	universal	129 8488 1081 000	

LumboTrain[°] | LumboTrain[°] Lady

DES	CRIPTION	SIZE	ITEM NO.
01	Nubbed pad, incl. cover	universal	129 8456 0071 000



LumboLoc[®] Forte

DESCRIPTION		SIZE	ITEM NO.
01	Sacral pad	universal	121 8404 0050 001
02	Lumbar pad	universal	121 8405 0050 001



SacroLoc[®]

DESCRIPTION		SIZE	ITEM NO.
01	Pads (2-part)	universal	119 8416 0000 000

DESCRIPTION		SIZE	ITEM NO.
01	CFRP frame	universal	129 8463 008#000
02	Velour pouch		1298461008#000
03	Velcro patches for frame		12984650000000
04	Abdominal pad (small)	universal	12985310000000



= Design (1) short, (2) long

SofTec[®] Lumbo

DESCRIPTION		SIZE°	ITEM NO.
01	Shell, short	1 – 3	129 8482 0031 00°
02	Shell, long	1-4	129 8482 0032 00°
03	CFR rods	universal	129 8483 007#000
04	Corset stay	universal	129 8484 008# 000
05	Abdominal pad (large)	universal	129 8481 0080 000
06	Elastic band		129 8485 008# 000
07	Inner knitted fabric, short	1 – 3	129 8486 0081 00°
08	Inner knitted fabric, long	1-4	129 8486 0082 00°
09	Outer knitted fabric, short	1 – 3	129 8487 0081 00°
10	Outer knitted fabric, long	1 – 4	129 8487 0082 00°



= Design (1) short, (2) long

DorsoTrain®

01 Reclinator, short universal 129 8456 0071 000 02 Reclinator long universal 129 8456 0072 000	DES	CRIPTION	SIZE	ITEM NO.
02 Reclinator long universal 129.8456.0072.000	01	Reclinator, short	universal	129 8456 0071 000
	02	Reclinator, long	universal	129 8456 0072 000



SofTec[®] Dorso

DES	DESCRIPTION		ITEM NO.
01	Reclinator (with pads)	universal	129 8498 1170 000
02	Reclinator pad	universal	129 8501 1000 000
03	Set A (cross + assembly kit)		129 8493 1000 000
04	Cross		129 8499 1170 000
05	Assembly kit		129 8502 1000 000
06	Set B (2 shoulder straps + 2 pads)	universal	129 8494 0080 000
07	Shoulder strap (contents: 2 straps)	universal	129 8496 0080 000
08	Shoulder pads (contents: 2 pads)	universal	129 8497 0080 000
09	Inner knitted fabric	1-4	129 8492 0080 00°
10	Outer knitted fabric	1-4	129 8491 008# 00°



ACCESSORIES & SPARE PARTS

= Design (1) short, (2) long

Spinova®

DESCR	IPTION	DESIGN	SIZE°
01	Abdominal fastening, set		
02	Base support with corset stays**		1 – 7
03	Fastening extension		
04	Tensioning straps, set**		1-5
05	Side clip, set (x10)		
06	Back setting device, set		
07	Shell system, set**	small	
07	Shelt System, Set	large	
08	Shell cover, set**	short	1-2
00	Shell Cover, Set	long	1-2
09	Bridging Frame	short	
07		long	
10	Prideing and	standard	
10	Bridging pad	short	
11	Cushion, bridging pad	short	
12	Core, bridging pad	short	
13	Sleeve, bridging frame, set (x4)		
14	Velcro patches for fastening the back setting device		
15	Velcro strips for fastening the frame, set (x10)		
16	Plastic rivets, bridging frame, set (x20)		
17	Shoulder pads, set		
18	Reclinator (without cover)	S-L	1-3
19	Immo screw set		
20	Immo Allen key		
21	Torx key for tensioning system		
22	Abdominal pad	large	
23	Torque wrench		



Stabi Classic	Support Plus	Support Plus Classic	Unload Classic	Unload Plus	Immo Plus	Immo Plus Classic	Immo	Osteo	ITEM NUMBERS
•	٠	•	•	•	٠	•	•		129 8545 0000 000
	•			•	•				129 8544 0000 00°
•	٠	•	•	•	٠	•	•		129 8565 0000 000
	•			•	•		•		129 8546 0000 00°
	٠			•	٠		•		129 8547 0000 000
	•			•	•		•		129 8548 0000 000
					٠	•	•		129 8554 0000 001
					٠	•	•		129 8554 0000 002
					•				129 8566 0001 00°
		•			•				129 8566 0002 00°
		•	•	•					129 8552 0001 000
		•	•	•					129 8552 0002 000
	•								129 8549 0000 000
	•								129 8549 0001 000
	•								129 8561 0001 000
	٠								129 8562 0001 000
			•	•					129 8553 0000 000
							•		129 8555 0000 000
			•	•		•			129 8550 0000 000
			•	•		•			129 8551 0000 000
								•	129 8517 0000 000
		•				•		•	129 8570 0000 00°
					•		•		129 8556 0000 000
					٠		•		129 8557 0000 000
	٠			•	٠		•		129 8558 0000 000
	٠			•	•				129 8481 0080 000
	٠			•	٠		•		129 8560 0000 000

° For patients taller than 170 cm

TREATMENT AREA: HAND, ELBOW AND SHOULDER

ManuTrain®

DES	SCRIPTION	SIZE	ITEM NO.
01	Velcro strip (for sizes 1–4)	1	129 5152 20*0 001
02	Velcro strip (for sizes 5–6)	2	129 5152 20*0 002
03	Stay, gray	universal	129 5151 1120 +00
- (Color (1) boigo (E) blue		

* = Color (1) beige, (5) blue + = Design (6) right, (7) left

ManuLoc[®] long

DESCRIPTION	SIZE	ITEM NO.	
01 Spare straps (x4), color: titan	universal	129 5343 0080 000	

ManuLoc[®] long Plus

DESCRIPTION	SIZE	ITEM NO.
⁰¹ Finger support with strap, color: tit	an 1	129 5342 0080 00#
⁰² Spare straps (x4), color: titan	universal	129 5343 0080 000

= Design (1) Size S, (3) Size L

ManuLoc[®] Rhizo

DES	CRIPTION	SIZE°	ITEM NO.
01	Thumb attachment	0-3	129 5346 0080 +0°
+ = [lesign (6) right, (7) left		



ManuLoc[®] Rhizo long

DES	CRIPTION	SIZE	ITEM NO.
01	Spare straps (x4), color: titan	universal	129 5343 0080 000

ManuLoc[®] Rhizo long Plus

DESCRIPTION		SIZE	ITEM NO.	
01	Finger support with strap, color: titan	1	129 5342 0080 00#	
02	Spare straps (x4), color: titan	universal	129 5343 0080 000	

= Design (1) Size S, (3) Size L









EpiTrain[®]

DES	CRIPTION	SIZE°	ITEM NO.
01	Strap, color: blue	0-6	119 6163 4080 00°

OmoTrain[®]

DESCRIPTION		ITEM NO.
01	Buckle, color: gray	119 7172 0000 000

OmoLoc[®]

DES	CRIPTION	SIZE	ITEM NO.
01	Shoulder strap	universal	129 7506 0000 000



01

SecuTec[®] Omo

DES	CRIPTION	SIZE	ITEM NO.
01	Cushion set and rivets	universal	129 7504 0000 000
02	Abdominal strap and fastening system	universal	129 7507 0000 000
03	Velcro strap set for arm support	universal	129 7508 0000 000
04	Small component set	universal	129 7509 0000 000
05	Hand support, cushion and rivets		129 7511 0000 +00
06	Shoulder strap	universal	129 7506 0000 000

+ = Design (6) right, (7) left





CUSTOM-MADE AND STANDARD PRODUCTS

MEASURING CORRECTLY FOR THE BEST TREATMENT SUCCESS

Medical aids can only be 100 % effective when they fit properly. This requires the exact measuring of the relevant body dimensions. That's why we have developed a range of order sheets for custom-made products. They help you to order the required Bauerfeind product quickly and easily in the design needed.

You can download the individual order sheets at any time using the following link: **www.bauerfeind.de/downloads**

We are happy to send you, free of charge, the "Train Active Supports" order pad for custom-made products (Order no. 19001060001002).



MEASURING ADVICE

Record all the measurements that have been specified on the measurement sheets – we can only ensure that the custom-made solution fits accurately if this data is complete.

Patients with a tendency toward pronounced swellings should be measured in the morning, if possible.

Indicate length markings on the skin using a kajal pencil or another skin-friendly pen to ensure reliable and reproducible circumference measurements.

If products are needed for both sides, please measure each side individually.

Do not apply tension when measuring.

After measuring, check whether all required lengths and circumferences have been recorded.

Knee

Measurements are taken with the patient standing, with a leg position of 30°, if the physical constitution allows.

Measurements should be taken when the leg is uncovered.

If the customer has weak connective tissue, skin flaps, etc., putting on pantyhose can make measuring easier.

Measure the side of the leg from the middle of the kneecap.

Depending on the product, the required minimum circumference dimensions must be noted.

For a thigh circumference of more than 60 cm, we recommend attaching a siliconized edge on the support.

Foot

Measurements are taken with the patient standing, if the physical constitution allows.

The foot length is measured on the outside, from the heel to the MTP joint of the little toe.

ORDER FOR CUSTOM-MADE PRODUCTS GenuTrain®



119

Identifier		
(e.g. patient name,		
abbreviation)		

	Circumfere	nce Length
	f ¹	e-f ¹
\vdash	f	e-f
$V \rightarrow V$	е	
	d	e – d
	c	e-c

All measurements must be entered to ensure proper treatment.

- $\begin{array}{ll} f^1 & \mbox{Thigh: Measure the circumference at least 14 cm above the middle of the knee (e). If the circumference at point f^1 is more than 47 cm, then the circumference must be measured at a higher point, where <math display="inline">e-f^1$ is at least 17 cm.
- f Base of thigh
- e Middle of knee
- d Narrowest point below the knee
- c Calf circumference: Measure the circumference at least 12 cm below the middle of the knee (e).

Notes:

- Take measurements with a slight flexion (30°).
- The length is measured on the outside of the leg. The distance between e and f¹ should not be shorter than the distance between e and c.
- If the circumference of f¹ exceeds 60 cm, we recommend attaching a silicone dotted band.

Fax +49(0)3662866-1699


 Quotation Order 	
Knee □ Right	🗌 Left
Color Titan	□ Beige
Extras	add 5 cm to the length)
□ Tightening strap*	
Velcro fastening*	□ Bottom
Patient details*** Male	🗌 Female
🗌 Obese	Pronounced muscles
Age	Height
*** Optional; please specify accurate fit of the finishe The purchaser is responsil	tion with a siliconized edge; not circular special anatomical features to ensure an
L Stamp	
Customer no.	
Signature	


Date

ORDER FOR CUSTOM-MADE PRODUCTS GenuTrain[®] P3



Identifier (e.g. patient name)





Notes:

- Take measurements while standing, with a slight flexion (30°).
- The lengths are measured on the outside of the leg.
- Please note:
- Min. circumf. f¹ = 39 cm
- Min. circumf. f = 37 cm
- Min. circumf. e = 30 cm
- Min. circumf. b¹ = 28 cm

Fax +49(0)3662866-1699

Quotation	🗌 Order
Knee □ Right	🗌 Left
Color Titan	□ Beige
Extras	dd 5 cm to the length)
Position of the correct Standard (lateral)	ive guide □ Alternative option (medial)
Patient details** Male Obese	 Female Pronounced muscles
Age:	Height: cm
length when measuring pat	ecial anatomical features to ensure an
The purchaser is respo size information is corr Custom-made solution	
All measurements mu	st be entered
to ensure proper treat	ment.
e Middle of knee	x. midpoint between e and f ¹) ow the knee, approx. 5–6 cm ference
Stamp	
Customer no.	
Signature	
Date	

ORDER FOR CUSTOM-MADE PRODUCTS GenuTrain[®] S



121

(e.g	entifier 1. patient name, previation)			Fax +49(0)366	28 66-1699
				Quotation	
		Circumference	Length	🗌 Order	
		f ¹	short: e – f ¹ = 18 cm long: e – f ¹ = 21 cm	Knee	□ Right
		f	e-f	□short	
					port, thigh $[e-f^1] = 18 \text{ cm}$)
	$[\cap]$	е		□ long (Length of supp	port, thigh [e-f¹] = 21 cm)
	$ \cup $			Extras	
	\searrow	d	e-d	Siliconized edg	ge* (add 5 cm to the length)
ſ				Patient details**	
		С	e-c	☐ Male	Female
		b1	e – b ¹ = 18 cm	Obese	Pronounced muscles
ł		0	e-b = 100111	Age	Height
All	Thigh		ensure proper treatment.	length when measu ** Optional; please spe accurate fit of the fir The purchaser is resp	cify special anatomical features to ensure an
			n above the middle of the knee (e).		
f	Base of thigh				
е	Middle of knee			I	I
d		nt below the knee			
С		ence: Measure the circ dle of the knee (e).	cumference at least 12 cm		
b ¹	Calf circumfer	ence 18 cm below the	knee		
T • T • T • T • T • f	Custom-made solutic Take measurements The length is measur The selection of the s Patient. The long vers	sion is recommended for pat um circumference measuren ee): at least 38 cm least 32.5 cm	color titan. on the individual measurements of the ients from a height of 1.80 m. nents are required for supplying the	Stamp Customer no.	
				Signature	
				Date	

ORDER FOR CUSTOM-MADE PRODUCTS GenuTrain[®] S Pro



(e.g.	ntifier patient nan eviation)	ne,		_ Fax +49(0)366	528 66-1699
				Quotation	
		Circumference	Length	🗌 Order	
		f1	short: e-f ¹ = 18 cm long: e-f ¹ = 21 cm	Knee - □ Left	□ Right
\vdash		f	e-f	_ □ short (Length of supp	port, thigh [e-f ¹] = 18 cm)
Y	Δ	e		□ long (Length of supp	port, thigh [e-f ¹] = 21 cm)
	\bigcirc	d	e-d	Extras	ge* (add 5 cm to the length)
		c	e-c	 Patient details** □ Male 	🗌 Female
		b ¹	e – b ¹ = 18 cm	- - Obese - Age	Pronounced muscles Height
f1	Thigh short version:	Measure the circumference 180	ensure proper treatment. orm above the middle of the knee (e). m above the middle of the knee (e).	length when measu ** Optional; please spe accurate fit of the fir The purchaser is resp	ecify special anatomical features to ensure an
f	Base of thi	igh			_
е	Middle of I	knee			
		point below the knee			
		nference: Measure the cir middle of the knee (e).	cumference at least 12 cm		
b1	Calf circun	nference 18 cm below the	e knee		
 Tak The The the The pro f (1) e (n) 	stom-made si ke measurem e length is me e selection of e patient. The e following m oduct: 4 cm above th middle of kne	long version is recommended fo		L Stamp Customer no. L	
C (. 2 cm betow t			Signature	
				Date	

ORDER FOR CUSTOM-MADE PRODUCTS MalleoTrain®



123

				Fax +49(0)3662	28 66-1699
	ntifier _ patient name,			Quotation	
	viation)			🗌 Order	
		Circumference	Length	Foot □ Right	🗌 Left
				Color	🗌 Beige
	_	b ¹	A-b ¹	Extras	
١		b	A-b	🗌 Velcro fastening	
				Measuring points	for pads**
		h		1	
				2	
		a	A–a (length of foot)	Patient details***	
			(10.1.9.11 01.1004)	🗌 Male	🗌 Female
-				🗌 Obese	Pronounced m
				Age	Height
				 Subject to additional Considered position of the 	l charges ne pads (position of ankle bi
				substantially from th	ne standard measurements
All m	easurements m	ust be entered to ensu	re proper treatment.	an accurate fit of the	cify special anatomical feat finished product.
b1	Base of calf: Mea above the heel (/	asure the circumference A).	at least 16 cm		onsible for ensuring that a Custom-made solutions c
b	Smallest circum	ference above the ankle	bone		
h	Instep (circular r	measurement of the inst	ep over the heel)		
а	Circumference b	ehind the MTP joint of th	ne little toe		
А	Heel				
A-a	Foot length with	out toe			

Notes:

- Only specify the pad measurement points if special positioning is required.
- Distance ${\boldsymbol \textcircled{O}}$ from the medial malleolus to the lateral malleolus (measured horizontally across the Achilles tendon)
- Distance ② from the medial malleolus to the lateral malleolus (measured over the sole of the foot)

Foot] Left				
Color Titan] Beige				
Extras	ening*					
Measuring po	oints for pa	ads**				
1						
2						
Patient detail] Female				
🗌 Obese] Pronoun	ced m	uscles	5	
Age		Height _				
 Subject to add Special positic substantially f Optional; pleas an accurate fit 	on of the pads rom the stand se specify spe	(position of a lard measure cial anatomic	ements)			ż
The purchaser is information is correturned.					ŧ	
Stamp						
Customer no.						
Signature						
Date						

ORDER FOR CUSTOM-MADE PRODUCTS EpiTrain[®]



Identifier (e.g. patient name abbreviation)	3, 		Fax +49(0)360			
/	Circumferenc	e Length	Arm □ Right	□ Left		
H	f ¹	e-f ¹	Color Titan 	☐ Beige		
$A \wedge$			Extras Siliconized edge* (add 3 cm to the length) 			
	f	e-f	Velcro fastening*	Bottom		
	e		Patient details** □ Male	∗ □ Female		
$\left[\right] $	d	e-d	Obese	Pronounced muscles		
	c	e-c	Age	_ Height		
Cull		I		nbination with a siliconized edge; not circular becify special anatomical features to ensure an		
All measureme	nts must be entered to e	ensure proper treatment.	information is correc	ponsible for ensuring that all size :t. Custom-made solutions cannot be		
f ¹ Upper arm: I crease of the		ce at least 10 cm above the	returned.			
f Base of the u	upper arm			_		
e Crease of the	e elbow					

- d Middle of the forearm muscle
- c Base of muscle: Measure the circumference at least 12 cm below the crease of the elbow (e).

Notes:

- Take measurements with a slight flexion (30°).
- The starting point to determine the length measurements is the crease of the elbow (e).
- The custom-made solution does not have a strap.

 Stamp				
Customer no.				
Signature				
Date				

QUALITY THAT CAN BE MEASURED – AND MADE TO MEASURE BAUERFEIND MEASUREMENT TECHNOLOGY



BODYTRONIC[®] 600 FULL-BODY SCANNER INCLUDING PRODUCT SUPPLY FOR THE LOWER EXTREMITIES

- Contact-free measurement for compression stockings, knee and foot supports as well as knee orthoses, with no need for special measurement clothing
- Optional upper body scan, export of 3D data
- Integrated leg volume measurement
- Standing space of 2.80 m²







BODYTRONIC[®] 610 FULL-BODY SCANNER INCLUDING PRODUCT SUPPLY FOR ALL EXTREMITIES

- Contact-free measurements for medical compression stockings, knee and foot supports, knee orthoses as well as for back, arm and shoulder products without the need for special measurement clothing
- Optional full-body scan, export of 3D data
- Integrated leg volume measurement
- Standing space of only 1.60 m²



SCANNER AND CAMERA Beams of light are projected to form patterns on the object being measured and are captured digitally. From the information in the images, the software generates a precise 3D model.



SELF-ROTATING MEASUREMENT PLATFORM Moves the patient steadily, thus excluding the potential for errors.



SELECTING THE PRODUCT WITH THE PERFECT FIT The BODYTRONIC 610 helps you to choose the perfect standard or custom-made product.

TERMS OF SALE AND DELIVERY OF BAUERFEIND AG (GENERAL BUSINESS TERMS) INTERNATIONAL

Prices

are quoted net ex works in Euro currency, plus freight and taxes, if any.

Dispatch

at buyer's cost and risk even if the buyer does not bear the transportation cost. Disruptions of operations affecting us or any of our suppliers, lack of suitable personnel, strike, lockout and force majeure are causes for extending delivery and relieve us of any delivery obligation for the time of the obstruction.

Returns

outside the warranty terms are only accepted after prior advise and approval of Bauerfeind to the return, provided the delivery documents are submitted; it is supposed that the returned goods are in proper state. We reserve the right of rejecting the acceptance of any returned consignment. We charge an administration fee of €10 for every consignment returned outside the warranty period. The administrational fee for goods returned 1 month or more after the date of delivery is additionally 20 percent of the net order value, subject to a maximum of €25 per article and piece; the fee is additionally 40 percent of the net order value of goods returned between 6 months and 12 months after the delivery date, subject to a maximum fee of €50 for each article and piece. Goods returned outside the warranty, especially custom made products, and 12 months or more after the date of delivery will not be accepted back.

Payment

to be made immediately after presentation of the invoice. First deliveries are only made against advance payment. If the customer is in default, any further delivery will only be made against advance payment. We charge interest on arrears at the rate of 8 percent above the applicable base rate from the time at which payment is delayed. Payments by check, draft, etc. are deemed to be made when the amount has finally been credited to our account.

Retention of title

We retain the title to all goods delivered until the buyer has paid all amounts outstanding under our business relations with him. The retention of title also applies proportionately to goods which have been processed or otherwise modified. If the goods are sold all claims on third parties are deemed to have been assigned to us without proviso up to the total amount owing to us. Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment we can demand that the buyer disclose the assigned claims and the names of the debtors, provides all information required for collection of the debt, delivers all required documents and informs the debtors (third party) of the assignment. The buyer is not permitted to sell, hypothecate or assign the goods as security. Hypothecation, confiscation or any other action of constraint taken by a third party shall be communicated to us by the buyer without

delay and he shall provide all information and documents to enable us to defend our rights. We undertake to release, on the buyer's demand, any security which exceeds a secured claim that has not been satisfied by more than 20 percent.

Complaints

Warranty claims due to obvious defects expire unless notified to us in writing without delay within the meaning of Article 377 HGB (German Commercial Code). For the rest, the time bar for warranty claims is 1 year after delivery of the goods purchased. If the goods at our choice. For reasons of hygiene, repairs can only be made to washed goods and goods which are of good hygienic condition. We will reserve the right of charging an adequate administration fee if unwashed goods are returned. If and to the extent to which we fulfil our duies of supplementary performance, in particular, the repair of defects or delivery of goods without defect, the customer cannot demand a price reduction or rescind the contract unless our attempts at supplementary performance have failed. Defects which are due to normal wear and tear, overstrain, misuse, neglect of care or non-compliance with the instructions for use are not covered by warranty.

Liability

for damage is excluded unless we or our servants are found to have acted with gross negligence or intention. This limitation of liability does not apply to damage due to injury of life, body or health or the violation of a major contractual duty.

Legal venue

The place of performance is the business seat of this company. The legal venue for all controversies from the business relations is Gera.

Partial invalidity

If any provision in these Terms of Sale and Delivery should be or become invalid or unenforceable, the other provisions will not be affected by this. The ineffective or unenforceable provision will be substituted by a valid and enforceable provision which reflects the original intention of the parties as closely as possible.

Data protection clause / call recording

For quality purposes, we randomly record calls between the buyer and our employees. The phone records exclusively serve the purpose of improving our performance and are deleted after evaluation. The data will not be disclosed to third parties. The buyer hereby consents to the records of telephone calls to that sole purpose. Consent is voluntary and may be revoked at any time in the future by the buyer.

Latest edition: February, 2016

A BETTER QUALITY OF LIFE FOR YOUR CUSTOMERS!

Bauerfeind AG

Triebeser Straße 16 07937 Zeulenroda-Triebes Germany P +49 (0) 36628 66-1660 F +49 (0) 36628 66-1699 E info@bauerfeind.com

AUSTRIA

Bauerfeind Ges.m.b.H. Hainburger Straße 33 1030 Wien P +43 (0) 800 4430-130 F +43 (0) 800 4430-131 E info@bauerfeind.at

BENELUX

Bauerfeind Benelux B.V. Waarderveldweg 1 2031 BK Haarlem

THE NETHERLANDS **P** +31 (0) 23 531-9427 **F** +31 (0) 23 532-1970 **E** info@bauerfeind.nl

BELGIUM

P +32(0) 2 527-4060 **F** +32(0) 2 792-5345 **E** info@bauerfeind.be

BOSNIA AND HERZEGOVINA

Bauerfeind d.o.o. Meše Selmovića 19 71000 Sarajevo P +387 (0) 33 710-100 F +387 (0) 33 619-422 E info@bauerfeind.ba

CROATIA

Bauerfeind d.o.o. Goleška 20 10020 Zagreb P +385(0) 1 6542-855 F +385(0) 1 6542-860 E info@bauerfeind.hr

FRANCE

Bauerfeind France S.A.R.L. B.P. 59258 95957 Roissy CDG Cedex **P** +33 (0) 1 4863-2896 **F** +33 (0) 1 4863-2963 **E** info@bauerfeind.fr

ITALY

Bauerfeind AG Filiale Italiana Piazza Don Enrico Mapelli 75 20099 Sesto San Giovanni (MI) P +39 02 8977 6310 F +39 02 8977 5900 E info@bauerfeind.it

NORDIC

Bauerfeind Nordic AB Storgatan 14 114 55 Stockholm P +46(0) 774 100 020 E info@bauerfeind.se

REPUBLIC OF NORTH MACEDONIA

Bauerfeind Dooel Skopje 50 Divizija 24 a 1000 Skopje P +389(0) 2 3179-002 F +389(0) 2 3179-004 E info@bauerfeind.mk

SERBIA

Bauerfeind d.o.o. 102 Omladinskih brigada 11070 Novi Beograd P +381 (0) 11 2287-050 F +381 (0) 11 2287-052 E info@bauerfeind.rs

SINGAPORE

Bauerfeind Singapore Pte Ltd. Blk 41 Cambridge Road #01-21 Singapore 210041 P +65 6396-3497 F +65 6295-5062 E info@bauerfeind.com.sg

SLOVENIA

Bauerfeind d.o.o. Dolenjska cesta 242 b 1000 Ljubljana P +386 (0) 1 4272-941 F +386 (0) 1 4272-951 E info@bauerfeind.si

SPAIN

C / San Vicente Mártir, nº 71-4º-7ª 46007 Valencia P +34 96 385-6633 F +34 96 385-6699 E info@bauerfeind.es

SWITZERLAND

Bauerfeind AG Vorderi Böde 5 5452 Oberrohrdorf P +41 (0) 56 485-8242 F +41 (0) 56 485-8259 E info@bauerfeind.ch

UNITED ARAB EMIRATES

Bauerfeind Middle East FZ-LLC Dubai Healthcare City Building 40, Office 510 Dubai P +9714 4335-684 F +9714 4370-344 E info@bauerfeind.ae

UNITED KINGDOM

Bauerfeind UK 85 Tottenham Court Road London W1T 4TQ P +44 (0) 121 446-5353 F +44 (0) 121 446-5454 E info@bauerfeind.co.uk

USA

Bauerfeind USA, Inc. 1230 Peachtree St NE Suite 1000 Atlanta, GA 30309 P +1800 423-3405 P +1404 201-7800 F +1404 201-7839 E info@bauerfeindusa.com

